

Determining the gravitational acceleration with a reversible pendulum

Objects of the experiments

- Measuring the oscillation periods T_1 and T_2 of a reversible pendulum for two suspension points.
- Tuning the reversible pendulum to the same oscillation period.
- Determining the gravitational acceleration from the oscillation period and the reduced length of pendulum.




Fig. 1 Oscillations of a reversible pendulum around the suspension points H_1 and H_2 .

Principles

Compound pendulum:

If a compound pendulum oscillates around its rest position with small deflections φ , the equation of motion is:

$$J \cdot \ddot{\varphi} + m \cdot s \cdot g \cdot \varphi = 0 \quad (I)$$

J : moment of inertia around the axis of oscillation,
 s : distance between the axis of oscillation and the centre of mass,
 g : gravitational acceleration, m : mass of the pendulum

The reduced length of the compound pendulum is defined as the quantity

$$s_r = \frac{J}{m \cdot s} \quad (II)$$

because its oscillation period

$$T = 2\pi \sqrt{\frac{s_r}{g}} \quad (III)$$

corresponds to that of a simple pendulum with the length s_r .

The moment of inertia J of the compound pendulum is, according to the parallel axis theorem,

$$J = J_S + m \cdot s^2 \quad (IV)$$

J_S : moment of inertia around the centre of mass axis

Therefore the reduced length of pendulum is

$$s_r = \frac{J_S}{m \cdot s} + s \quad (V)$$

Reversible pendulum:

The reversible pendulum is a particular type of the compound pendulum. There are two edges H_1 and H_2 that allow to choose the suspension point. Two masses $m_1 = 1000$ g and $m_2 = 1400$ g on the straight line H_1H_2 can be shifted so that the oscillation period is tunable. The goal of the tuning is to achieve equal oscillation periods around both edges. In this case, the reduced length of pendulum is equal to the distance $d = 99.4$ cm between the edges. This latter statement can be understood from the following consideration:

- Hang the reversible pendulum on the edge bearing with the edge H₂, and measure the period 50 · T₂.
- Slide the mass m₂ to the position x₂ = 55 cm, and measure 50 · T₂ at first and then 50 · T₁.
- Slide the mass m₂ towards the edge H₂ in steps of 5 cm; each time measure the two oscillation periods. Plot T₁² and T₂² as functions of x₂ and, if necessary, repeat the measurement of the oscillation periods.
- Next slide the mass m₂ towards the edge H₁ in steps of 5 cm starting from x₂ = 45 cm, and measure the two oscillation periods each time.

Evaluation and results

The two measured curves T₁²(x₂) and T₂²(x₂) intersect near x₂ = 30 cm and x₂ = 65 cm (see Fig. 4). The enlarged sections in Figs. 5 and 6 shown that the curves intersect at T² = 4.039 s² and T² = 4.014 s², respectively. With the mean value T² = 4.027 s² Eq. (X) gives

$$g = \frac{4 \cdot \pi^2 \cdot d}{T^2} = 9.74 \frac{\text{m}}{\text{s}^2}$$

Measuring example

Table 1: Oscillation periods T₁ and T₂ around the edges H₁ and H₂, respectively, as functions of the distance x₂ between the mass m₂ and the edge H₁.

x ₂ cm	50 · T ₁ s	T ₁ ² s ²	50 · T ₂ s	T ₂ ² s ²
20	106.0	4.494	101.9	4.153
25	103.1	4.252	101.2	4.097
30	100.8	4.064	100.6	4.048
35	99.4	3.952	100.1	4.008
40	98.8	3.905	99.8	3.984
45	98.2	3.857	99.6	3.968
50	97.9	3.834	99.5	3.960
55	98.3	3.865	99.6	3.968
60	99.1	3.928	99.8	3.984
65	99.9	3.992	100.0	4.000
70	101.1	4.088	100.7	4.056
75	102.2	4.178	101.7	4.137
80	103.2	4.260	102.2	4.178
85	104.8	4.393	103.6	4.293




Fig. 5 Enlarged section of Fig. 4 around x₂ = 30 cm with a non-linear interpolation of the measured values. Point of intersection: (x₂ = 31 cm, T² = 4.039 s²).

Fig. 4 Squared oscillation periods around the edges H₁ (●) and H₂ (■) as functions of the distance x₂ between the mass m₂ and the edge H₁.




Fig. 6 Enlarged section of Fig. 4 around x₂ = 65 cm with a non-linear interpolation of the measured values. Point of intersection: (x₂ = 66 cm, T² = 4.014 s²).

