
Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

8/95-kem/Sf-

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**Demonstration-experiment frame**  
**Shelf**  
**Profile rail**  
**Console**

Fig. 1

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Demonstration-experiment frame, mounted on large trolley (598 11) by means of two bench clamps with pin, 301 05 (598 11 and 2x 301 05 not included in scope of supply)

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The demonstration-experiment frame and supplementary equipment allow assembly of clear, comprehensible experiment setups in which all apparatus are visible and e.g. the wiring of electrical circuits is easy to follow. Power supplies, experiment apparatus and measuring instruments can be arranged so as to be observable by all students. Display and measuring instruments in particular can be placed or suspended at a great height above the surface of the laboratory bench.

When the demonstration-experiment frame is mounted on a trolley, experiments can be prepared in advance and transported to the classroom.

## 1 Safety notes

Only attach the profile rails ⑨ in the frame so that the deep slot with the brush strips (see enlarged detail in Fig. 2) is facing downward, to ensure that the apparatus can be mounted securely.

When using on a trolley, attach the frame securely to the two holes ⑩ in the trolley surface using two bench clamps with pin (301 05).

Be sure that there is enough headroom in door frames for your experiment setups.

Avoid abrupt movements of the trolley to which an experiment setup is mounted. If this cannot be assured (e.g. due to knobbed floor surfaces or other irregularities), make sure that none of the apparatus can fall out of the assembly.

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Exploded drawing of demonstration-experiment frame (301 300)  
assembly of parts in the order of their numbers recommended  
Inset drawings: three different possible assemblies

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### 301 300 Demonstration-experiment frame (Fig.2)

- ① 1 Transverse bar, l = 93 cm
- ② 3 Thumb screws, for mounting support rods in the transverse bar holes
- ③ 2 Side sections, 35 cm x 91 cm  
For attaching the side sections to the transverse bar:
- ④ 4 Spring washers
- ⑤ 4 Allan (hex) screws M5 x 25
- ⑥ 1 Allan key
- ⑦ 8 Covers  
For mounting the profile rails:
- ⑧ 4 Knurled screws
- ⑨ 2 Profile rails, 93 cm x 4.5 cm x 2 cm
- ⑩ Holes for mounting the frame to a trolley using bench clamps with pin (301 05, not included in scope of supply)

Dimensions (assembled):

Height: 91 cm  
Width: 97 cm  
Depth: 35 cm

### 301 311 Profile rails

Supplementary equipment for demonstration-experiment frame  
Dimensions: 93 cm x 4.5 cm x 2 cm  
with two knurled screws

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Abmessungen: 30 cm x 16 cm

### 301 310 Shelf

Supplementary support surface for demonstration-experiment  
frame. The shelf can be inserted in demonstration-experiment  
frames either without profile rails or having one or two profile  
rails.

Dimensions: 30 cm x 97 cm x 2.5 cm

With four knurled screws for sturdy mounting (only necessary  
when no profile rails are used.)

### 301 312 Console

For attachment to a profile rail as support surface in front of the  
frame, e.g. for an oscilloscope.

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## 3 Assembly

Assemble the demonstration-experiment frame as shown in  
Fig. 1. Using the Allan key, make sure that all Allan screws are  
screwed tight. In particular, make sure that the profile rails are  
attached with the groove facing upwards (slot with brush strips  
facing downwards).

When using with a trolley (e.g. 598 11), fasten the demonstra-  
tion-experiment frame securely using two bench clamps with pin  
(301 05).

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
Fig. 3

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#### Demonstration-experiment frame without supplementary equipment

In this setup, the distance from the benchtop to the lower  
profile rail is 22.5 cm, providing sufficient space for supply  
units in large housings

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#### Demonstration-experiment frame with adhesive magnetic board

In this setup, the distance from the benchtop to the lower  
profile rail is 9 cm, providing sufficient space for storage trays

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Demonstration-experiment frame with one shelf.  
The coil is mounted in the holder for tubes and coils (301 340)

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Demonstration-experiment frame with one profile rail and three shelves

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Demonstration-experiment frame with three profile rails, three shelves, one console and adhesive magnetic board.  
In the configuration shown here, the adhesive magnetic board serves as a neutral background; however, it can be used immediately just by turning the trolley.  
The plug-in board is attached using the plug-in board carrier.

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Laminated 1 Hz oscillator circuit with coil with high inductivity (517 011) and capacitor 40  $\mu$ F (517 021)

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Refraction at a lens  
with apparatus from the equipment sets 463 701 and 463 71

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Forces at pulleys  
with apparatus from the equipment set 343 001

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### 4 Additional accessories (examples)

- 301 301 Adhesive magnetic board  
For magnetic attachment of apparatus, e.g. from  
optics and mechanics.  
Dimensions: 93 cm x 62 cm
- 301 320 Plug-in board carrier  
For attaching plug-in boards A4 or A3
- 301 340 Holder for tubes and coils  
For positioning round apparatus in front of the  
demonstration-experiment frame  
Dimensions: 40 cm x 29.7 cm  
Midpoint spacing of apparatus: 14 cm
- 666 470 CPS-holder with bosshead, height-adjustable  
For attachment of apparatus with stand rod or by  
means of a universal clamp.