Fig. 1

Wurfgerät
Projection apparatus

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The apparatus is intended for investigation of the laws of motion for vertical, oblique and horizontal projection, and for demonstration of the mutual independence of the motions.

Literature:
Central card file, mechanics (598 611).


1 Description

The projector consists of a metal sleeve with a rectangularly attached metal rod of approx. 12 mm diameter for setting it up on stands. The acceleration rod is arranged axially in the sleeve. It can be pressed back against the force of a strong helical spring and retained in this position by means of a thumb screw. Two plastic balls each with one hole are used as projectiles. They can be fitted on the ends of the tensioned acceleration rod.

2 Operation

Set up the projection apparatus either according to fig. 2 (horizontal projection, independence principle) or to fig. 3 (oblique projection) in a stand. Tighten all screws.

Fig. 2

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To set reproducible spring forces, use the boss (301 01) as a defined stop for the end of the acceleration rod (see fig. 3).


Fig. 3


Fig. 4

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Press acceleration rod (a) (see fig. 4) in the arrow direction to tension the spring when the screw (b) is undone, and then tighten screw (b); (greatest spring tension imparting the maximum initial speed, with acceleration rod (a) moved up to the stop).

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After loading the springs, place the projectiles on the end (c) (or on both ends to demonstrate the independence principle) of the acceleration rod; quickly undo screw (b).