
6/96-Wiß/Sf-

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These pieces of equipment are intended for dynamics experiments on the precision metal rail (460 81 or 82). The motion of the trolley over time is recorded e.g. on metallized paper with the time recorder (337 18).

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1 Safety note

Do not vary the screws ④ for adjusting the bearing friction in a careless or random manner (overtightening will damage the mechanism, and loosening completely dismounts the wheel axles).

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2 Scope of supply, description, technical data

- ① to ④ Trolley 1 (337 00)
with 4 wheels on point bearings, which only touch the surface of the precision metal rail (460 81).
Weight: $85 \text{ g} \pm 1 \text{ g}$
Dimensions: $10.6 \text{ cm} \times 5.6 \text{ cm} \times 4.2 \text{ cm}$
- ① 4-mm sockets (one each on the front and rear) for attaching of spring and buffer (337 03)
- ② Plug-in holder for slotted weights ⑥ (from 337 04) and additional weight ⑦ (337 05)
- ③ Holding clip for metallized recording paper (337 39 or from 337 18)
- ④ Bearing screws (on both wheel axles) for readjustment of bearing friction (see section 4).
- ⑤, ⑥ Set of driving weights (337 04)
- ⑤ Flat weight, $5.2 \text{ g} \pm 0.1 \text{ g}$ for accelerating the trolley with a force of 0.05 N (see Fig. 2)

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- ⑥ Slotted weight (3x), $5\text{ g} \pm 0.1\text{ g}$, for attachment on to ② to increase the weight of the trolley to 90 g, 95 g or 100 g and for use in connection with ⑤ to increase the accelerating force to 0.1 N, 0.15 N and 0.2 N.
- ⑦ Additional weight, 100 g, for attachment on to ②

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3 Operation

Additionally required equipment:

- 1 Precision metal rail, 1 m long 460 81
- 1 Time recorder with metallized recording paper 337 18

For acceleration with driving weights:

- 1 Single deflection pulley on clamp rider (with Perlon thread) 337 14

For collision experiments and experiments on momentum

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Setup for investigation of uniformly accelerated motion $s = \frac{1}{2} a t^2$ and for derivation of Newton's equation of motion $F = m \cdot a$.


Fig. 3

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$$m = m_{\text{Wagen}} + m_{\text{Antrieb}}$$

Note: When using the driving weights ⑤ and ⑥ (337 04), the slotted weights ⑥ should be distributed according to Fig. 3 between trolley and weight ⑤ so that the total accelerated mass m remains constant:

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4 Maintenance

Maintenance-free bearings! Do not oil!

Readjustment of bearing friction is only required in those cases where play of the wheels has changed due to external stresses. In order to achieve an optimum friction value, loosen screws ④ by means of a small flat pliers until just a little play of the axles is noticeable (see section 1).