
5/95-d.R-

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The metal wheel and axle is used to produce and compare torques and to prove that the point of application of a force can be displaced along the line of action of force. It can also be used to demonstrate the mode of operation of a belt drive with different transmission ratios.

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1 Scope of Delivery; Technical Data

- ① Metal disk with three steps
Dia. of steps: 5 cm, 15 cm, 20 cm (1:3:4)
Weight: 0.7 kg
- ② Wheel axle with ball bearing and knurled-head screw which can be inserted from both sides
- ③ 4 rows of boreholes (distance from centre and spacing 4 cm) for inserting the hollow pins ⑦ (points of application of load or force respectively)
- ④ 4 boreholes to insert the crank ⑥
- ⑤ Grooves for inserting cords
- ⑥ Crank which can be inserted in boreholes ④
- ⑦ 4 Klemmstifte matching boreholes ③
- ⑧ Fishing line (309 48)

2 Operation

Additionally required equipment:

1 Dynamometer, 5 N	314 161
10 Weights, 50 g,	from 342 61
Stand material	
For demonstrating the belt drive:	
1 Second metal wheel and axle	342 75
or	
1 Pulley with hooks and rod	341 65

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Important:

Wind the fishing line once or twice around the step used for the measurement before suspending the weights or the dynamometer respectively.
Make the knots for the fishing line loops *small* and do not place them into the grooves.
Prior to starting the experiment, set the dynamometer in position of use to zero.

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When setting up the apparatus as *dissymmetrical lever* (Fig. 2), make sure that the points of application A and B of load and force are at the same level.

For *producing and comparing torques* (Fig. 3), fasten the fishing lines on the pins of two wheel steps with the radii r_1 and r_2 , insert the cords into the respective grooves and suspend weights.

For *lifting a load* (Fig. 4), place a cord with suspended weight into one of the grooves and turn the crank. For quantitative experiments a dynamometer can be used for pulling at the crank.

To demonstrate that the *point of application of force can be displaced along the line of action of the force* (Fig. 5), fasten the cord with suspended weights or dynamometer on the hollow pins and insert these into the respective boreholes. It is advisable to use only the boreholes up to the height of the axle or underneath to avoid an unstable position of equilibrium.

To demonstrate the *belt drive* (Fig. 6), insert the crank handle from the plain side. Transmission ratios are demonstrated by placing a loop of a cord successively into the grooves of the different steps.