08/90

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The instrument is used for demonstrating capillarity effects.

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1 Safety Notes

- Do not fill in acids!
- Be careful during experimentation with mercury! Carry out the experiment using the large plastic tray (309 63). Remove any spilt mercury using mercury adsorbent (306 83).

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2 Description

Two side walls of the wedge-shaped glass trough are sized about 105 mm x 40 mm, the third side about 15 mm x 40 mm. The wedge angle amounts to approx. 7°.

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3 Experiments

Before each experiment the trough must be well cleaned and dried. Then it is first filled with a non-wetting liquid, e.g. mercury. The mercury surface will be lower in the acute angle of the wedge than at its base; the smaller the distance between the two plates, the lower will be the mercury surface (capillary depression).

On the other hand, a liquid having a low surface tension, such as water, will rise upwards in the narrow space between the plates (capillary rise). The narrower the space between the plates, the higher the liquid will rise.

If the horizontal direction is called the x-direction and the vertical direction is the y-direction, then we have the following approximate relation for the co-ordinates (x, y) of any arbitrary point of the curve which the surface of the liquid describes:

$$x y = \frac{\sigma}{\rho g \cdot \operatorname{tg} \frac{\beta}{2}}$$

where

σ = coefficient of surface tension

ρ = density of the liquid

g = acceleration due to gravity

β = angle formed by the two glass plates

Now the expression on the right-hand side of the equation will have a certain value for each liquid, which will be constant for this particular liquid. Therefore we have

$x y = \text{const.}$

This is the equation of an equilateral hyperbola, of a curve which limits the surface of the liquid. The shape of the hyperbola will depend on σ and ρ of the liquid, as seen from the first equation. If, for example, the surface tension of water is lowered by addition of a wetting agent (such as is supplied, for example, with the wave trough, 408 02), then a change in the shape of the curve is distinctly observed.

Fig. 1


Fig. 2