
10/95-Sf-

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Verstellbarer Spalt Adjustable Slit

Fig. 1

The adjustable slit is used in all optical experiments where a continuously and reproducibly variable slit is required, particularly for diffraction and interference experiments and for spectral analysis.

1 Safety note

Impurities (e.g. grains of dust) adhering to the slit edges can cause undesirable diffraction effects; do not remove these using a metal object; use moist, lint-free paper or a fine moistened brush.

Always loosen screw ④ before realigning the slit.

2 Description, technical data

- ① Slit with chamfered edges:
Width: 0 to approx. 2.5 mm, continually adjustable
Height: 20 mm
- ② Lever for adjusting the slit width
- ③ Slit holder, rotatable by 360° in frame ⑥
- ④ Arresting screw for slit holder ③
- ⑤ Scale with index mark to indicate the slit width
Graduation: 0 to 1.6 mm with 0.2 mm divisions
Accuracy: 0,02 mm
- ⑥ Frame (dia. 13 cm)
on rod (8 cm x 1 cm diam.)

3 Operation

Set up the apparatus in the experiment arrangement so that the bevelled edges of the slit face the light source and the light beam is perpendicularly incident on the slit.

To optimally utilize the light, align the slit and the lamp filament parallel to each other. If necessary, loosen screw ④ and turn holder ③; remember to retighten the screw.

Rough adjustment of the experiment setup is best carried out when the slit is opened wide. Set the slit to the desired width using lever ② before fine adjustment of the setup.