
9/95-kem-

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Fig. 1

This apparatus can be used to demonstrate refraction, reflection and total internal reflection of light at the water-air boundary at various angles.

1 Safety notes

Before removing the diaphragm insert ③ (e.g. when replacing the halogen lamp), allow it to cool; after longer periods of operation, the temperature inside can reach over 100 °C!

Do not supply the halogen lamp with voltages greater than 12 V.

When using fluorescein sodium (469 74), rinse out the transparent plastic water cell with pure water afterwards.

Carefully wipe the water cell dry after each use to prevent scale build-up.

Do not clean the apparatus using aggressive fluids such as e.g. solvents or alcohol. We recommend using water with just a little washing-up liquid added.

Do not touch the bulb of the halogen lamp with your fingers. If you should accidentally leave fingerprints on the lamp, remove these immediately using a lint-free rag moistened with alcohol to prevent the fingerprints from being burned into the quartz glass.

2 Description, technical data

- ① 4-mm connector sockets for supplying the halogen lamp ②
- ② Halogen lamp: 12 V/20 W (309 04 794)
- ③ Rotating diaphragm insert with two slit diaphragms: 7 cm dia.
Number of slits: 1
10; angle between light beams: 10°
- ④ Transparent plastic water cell:
Dimensions: 24 cm x 13 cm x 2,5 cm
Filling volume: approx. 700 ml
- ⑤ Water-level mark
- ⑥ Drain outlet with silicone hose and hose clamp
- ⑦ Panel: 30 cm x 29.7 cm

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3 Operation

Recommended voltage source: e.g. transformer 6/12 V; 30 W (562 73)

Place the apparatus in a demonstration experiment frame (301 300) between two profile rails or mount it with the pair of stand feet (301 339).

Carry out the experiments in a darkened room.

Fill the cell with water approximately up to the mark ⑤ and connect the halogen lamp.

When the diaphragm is configured so that only one light beam emerges, you can continuously vary the angle at which it meets the water-air boundary by turning the diaphragm insert.

The 10-fold slit diaphragm allows you to observe ten light beams which are incident on the water-air boundary simultaneously.

4 Replacing the halogen lamp (see Fig. 2)

Caution! First allow the diaphragm insert ③ and lamp ② to cool!

Demount the diaphragm insert (see Fig. 2) and remove the defective lamp.

Place the new halogen lamp (459 01), still in the protective envelope, in the socket and carefully turn it until it stops. Only then should you remove the protective envelope.

Remount the diaphragm insert.


Fig. 2