5/90

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1 Description

The copy of a Rowland grating is copied from an original Rowland grating on a transparent collodium foil and is about 23 mm x 35 mm in size. The grating has about 600 lines per mm; the grating constant, the distance between two adjacent lines, is therefore 1/600 mm. The foil is inserted in 50 mm x 50 mm metal frame between two glass plates.

The collodium foil is subject to minor changes of shape due to meteorological influences. Therefore, the spectroscopical resolution is not quite the same as in case of the original grating.

2 Use

The frame is clamped into the holder with spring clips (460 22) and either mounted on a saddle base (300 11) or on an optical bench (e.g. 460 43).


3 Experiments

a) Projection of diffraction spectra
In order to produce brilliant diffraction spectra for group viewing we recommend using the mercury vapour lamp (451 15) in an experimental arrangement of the type shown in Fig. 2. The lens, $f = 200$ mm, forms an image of the slit on a wall 2 m to 3 m away. After mounting the grating directly behind the projecting lens the spectral lines of mercury can be seen on the wall in a bright display.

Fig. 2

- ① 451 15/19
- ⑤ 460 22
- ⑩ 460 04
- ⑪ 460 43
- ⑬ 460 14
- ⑭ 471 23
- ⑰ 460 02
- ⑲ 451 30

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b) Individual viewing of spectra using spectrum lamps

The spectra of spectrum lamps (451 011 - 451 111 with 451 16 and universal choke 451 30, see relevant Instruction Sheet) are best viewed in an arrangement of the type shown in Fig. 3. If the eye is brought up close to the grating, spectra of several orders can be viewed towards both sides.


Fig. 3

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c) Measurements of wavelengths or of grating constants

These measurements are e.g. done in an experimental arrangement as shown in Fig. 4, with a spectral tube which supplies a line spectrum with defined wave lengths. Please see to it that the grating is aligned in parallel to the rule at a distance of about 70 cm to 120 cm from it. The spectral tube must be directly in front of the rule.

The angle of diffraction follows from the distance between the rule and the grating and from the deviation of the spectral lines that can be read off on the rule. If the wavelengths are known, the grating constant must be computed and if the grating constant is known, the wavelengths can be determined.


Fig. 4