10/1986

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Satz spannungsoptischer Modelle Set of Photoelastic Models

The photoelastic models serve to demonstrate the double refraction which occurs when mechanical stresses are produced as a result of the elastic deformation of transparent solids.

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1 Safety Notes

- Protect the models against heat (for example direct sun light); max. permissible temperature: 90 °C.
- Do not subject the models to continuous loads.
- Do not violently bend the models.
- For cleaning use only non-aggressive cleaning agents.

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2 Description

The set consists of a simple straight rod, a hook, a right angle, and a right angle with its vertex drilled out.

The transparent models are made of unsaturated synthetic polyester resins whose photoelastic sensitivity is high. They have been stress-relieved by a suitable heat treatment and maintain this condition over long periods. Under the influence of external forces, stresses are set up in them which affect polarized light and under suitable conditions become visible as coloured fringes.

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3 Principle

These fringes are due to mechanical stresses causing anisotropy in the previously isotropic solid. It then exhibits double refraction and splits the incident polarized light into two components whose planes of oscillation are perpendicular to one another. Their phase velocities will differ by an amount increasing with the stresses in the body. The emergent light is therefore elliptically polarized and the analyzer can no longer extinguish it everywhere. As double refraction depends very strongly upon the wavelength, the system of fringes will be coloured if white light is used. The colours are mixed colours comprising those parts of the spectrum which are not extinguished at any given part of a fringe.

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Observation of the models in plane-polarized light yields a rather complicated picture because, apart from the coloured fringes due to the stresses in the material, there will be additional dark regions which depend upon the position of the model in the path of the rays. If the principal stresses of the body are parallel to the planes of oscillation of the polarizing filters, then double refraction cannot become effective and the field will remain dark even when the models are stressed. The body must, therefore, be rotated relative to the planes of polarization of the two filters.

The dependence of the pattern obtained upon the position of the body is avoided, and the pattern is more easily understood, if circularly-polarized light is used, and, therefore, this kind of light is preferred for the observation of photoelastic effects. Illuminated and coloured zones then correspond unambiguously, and independently of the direction of stress, to particular states of stress.

3 Assembly

3.1 Simple arrangement

These stresses can be observed qualitatively if the models are simply held between two crossed polarizing filters and bending them is attempted.

3.2 Arrangement using circularly polarized light

The items of apparatus are listed in the order in which they are incorporated into the assembly (see fig. 2):


Fig. 2

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Now the polarizing filters are rotated to 45° and 135°. If this results in a bright field of view, one of the $\lambda/4$ plates should be rotated through 90°.

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5 Experiments

The models are held in the parallel part of the beam.

The plain rod is held between both hands and an attempt is made to bend it in a plane which is perpendicular to the direction of the beam. The black, neutral thread and the parallel, coloured lines of equal stress, the isochromatics, are clearly seen; the latter provide a method of calibrating the process.

The limbs of the plain right angle are compressed or pulled apart by both hands, when the stress lines will converge to a singular point in the vertex.

The right angle with the hole will not exhibit this singularity when subjected to similar stress.

The hook can be pulled strongly with one finger or clamped to a stand by means of a Leybold multiclamp and loaded with a weight.