

Forces and work  
Force-transforming apparatus

## Forces and paths at a fixed pulley - Stand setup

## Objects of the experiment

1. Investigating the relation between load and force
2. Investigating the relation between the load path and the force path

## Setup



## Apparatus

1 Pulley, 100 mm diam., plug-in, pair, set of 2	340 921ET2
1 Weight, 0.5 kg	315 38
1 Weight, 0.2 kg	683 11
1 Weight, 0.1 kg	683 10
1 Precision dynamometer, 5 N	314 161
1 Metal rule, 0.5 m	460 97
2 Pointer, pair	301 29
3 Stand bases MF	301 21
1 Stand tube, 1000 mm, 10 mm diam.	608 052
1 Stand tube 750 mm, 10 mm diam.	608 051
1 Stand tube 450 mm, 10 mm diam., set of 2	666 609ET2
1 Stand rod, 250 mm, 10 mm diam.	301 26
2 Support block	301 25
1 Support clip, for plugging in	314 04ET5
1 Demonstration cord	309 50

## Carrying out the experiment

1. Relation between load and force:
  - Using the dynamometer, measure the gravitational forces of the three weights one after another, and enter them into the table as load  $F_1$ .
  - Then assemble the experimental setup as shown in the illustration.
  - Make the zero adjustment at the unloaded dynamometer.
  - Suspend the weight 0.5 kg.
  - Read the force  $F_2$  at the dynamometer.
  - Then suspend the other two weights one after another and read the corresponding forces.
2. Relation between the load path and the force path:
  - Set the load path  $s_1$  with the pointers at the left stand rod to a fixed value of 10 cm.
  - Align the load so that its lower edge is at the height of the lower pointer.
  - Mark the position of the dynamometer with the upper pointer at the stand rod in the middle.
  - Loosen the screw at the support block and lower the dynamometer until the load has passed the fixed path  $s_1$ .
  - Fasten the screw at the support block.
  - Mark the new position of the dynamometer at the stand rod in the middle with the lower pointer.
  - Measure the distance between the upper and the lower pointer at the stand rod in the middle, and enter it into the table as force path  $s_2$ .
  - Repeat the experiment with other load paths.

## Measuring example

Body	*Load $F_1$ in N	*Force $F_2$ in N
weight 0.5 kg	5	5
weight 0.2 kg	2	2
weight 0.1 kg	1	1

\*rounded values

Load path $s_1$ in cm	Force path $s_2$ in cm
10	10
15	15
20	20

## Evaluation

1. At the fixed pulley, the load and the force are equal:  $F_1 = F_2$ . Only the direction of the force is changed.
2. The load path and the force path are equal:  $s_1 = s_2$ .

Remark:

The measured values can be used to calculate the mechanical work performed at the fixed pulley:

$$W = F_1 \cdot s_1 = F_2 \cdot s_2$$