

Motions


Inertia and Newton's law

Relationship between acceleration, force and mass -  
Measurement using an electronic stopwatch

Objects of the experiment

1. Measuring the time  $t$  required by a body of mass  $m_1$  or  $m_2$ , respectively, to cover a certain path  $s$  if the force  $F$  acting on the body is changed
2. Calculating the acceleration  $a$  of the body
3. Representing the relation between acceleration and force in an  $a$ - $F$ -diagram

Setup



Apparatus

1 Track, 1.5 m .....	337 130
1 Trolley .....	337 110
1 Additional weights, pair .....	337 114
1 Holding magnet .....	683 41
1 Holder for combination spoked wheel .....	337 463
1 Combination spoked wheel .....	337 464
1 Combination light barrier.....	337 462
1 Slotted mass hanger, 10 g, small.....	315 410
4 Slotted weights, 10 g, red .....	315 416
1 Electronic stop-clock P .....	313 033
1 Connecting leads, 19 A, 50 cm, black, pair .....	501 451
1 Multi-core cable, 6-pole, 1.5 m.....	501 16
1 Fishing line, set of 2.....	309 48ET2
1 Plug-in axles, set of 2 .....	340 811ET2

- Repeat the experiment.

- Calculate the accelerations  $a$  from the quotients  $\frac{2s}{t^2}$ .


Measuring example

\* $m_1 = 1050$  g, \* $m_2 = 1550$  g

Path $s$ in cm	50	50	50	50
*Force $F$ in N	0.2	0.3	0.4	0.5
Time $t_{m_1}$ in s	2.53	2.02	1.72	1.53
Time $t_{m_2}$ in s	3.01	2.48	2.10	1.87
Acceleration $a_{m_1}$ in $\text{cm/s}^2$	15.6	24.3	34.4	43.5
Acceleration $a_{m_2}$ in $\text{cm/s}^2$	11.1	16.1	22.7	28.6

\*Force  $F$ : round values

Evaluation



In the case of uniformly accelerated motion, the acceleration varies proportionally to the acting force:  $a \sim F$ .

The greater the mass of a body, the smaller its acceleration if the force is the same.

Carrying out the experiment

- Adjust the voltage at the holding magnet so that the trolley with the additional weight is just held.
- Define the starting point with the movable interrupter flag on the trolley, and read it from the scale of the track.
- Position the light barrier at a distance of 50 cm from the starting point.
- Release the motion by pressing the START/STOP key at the stopclock.
- Wait until the interrupter flag passes the light barrier, and read the time from the stopclock.
- Reset the stopclock to zero by pressing the RESET key.
- Enhance the accelerating force step by step by putting the slotted weights 10 g from the trolley on slotted mass hanger one after another.
- Repeat the measurement for each force.
- Enhance the mass of the trolley by putting another additional weight on it.

