


## Determining the wind speed with a Prandtl pressure probe – Measuring the pressure with a pressure sensor and CASSY

### Objects of the experiment

- To measure the dynamic pressure.
- To determine the wind speed.
- To record pressure profiles from different distances.

Fig. 1: Prandtl pressure probe for measuring the static pressure  $p_s$  and the total pressure  $p_{tot}$ .



### Principles

The Prandtl pressure probe used in this experiment allows measuring the total pressure  $p_{tot}$  (head opening positioned against the direction of flow) and the static pressure  $p_s$  (head opening positioned across the direction of flow).

The wind speed  $v$  can be determined by measuring the pressure difference  $\Delta p$  while the Prandtl pressure probe is pointing against the direction of flow. In this case the pressure difference is identical with the dynamic pressure  $p_d$ :

$$\Delta p = p_d = p_{tot} - p_s \quad (I)$$

The average wind speed can now be calculated with the following equation:

$$v = \sqrt{\frac{2}{\rho} \cdot p_d} \quad (II)$$

Density of the air:  $\rho = 1.2 \frac{\text{kg}}{\text{m}^3}$

*Remark: In this experiment the Prandtl pressure probe is used in combination with CASSY and the Pressure sensor S,  $\pm 70$  hPa (Fig. 2). A precision manometer is used in P1.8.5.3.*

**Apparatus**

1 Suction and pressure fan .....	373 041
1 Prandtl pressure probe.....	373 13
1 Pressure sensor S, $\pm 70$ hPa .....	524 066
1 Sensor-CASSY 2 .....	524 013
or	
1 Mobile-CASSY .....	524 009A
or	
1 Pocket-CASSY 2 Bluetooth.....	524 018
1 CASSY Lab 2 .....	524 220

*Additionally required: 1 PC with Windows XP or higher*

*Optional:*

1 Stand base, V-shaped, small .....	300 02
1 Stand rod, 25 cm, 12 mm $\varnothing$ .....	300 41
1 Leybold multiclamp .....	301 01

**Setup**

- Equip the suction and pressure fan with the small nozzle ( $\varnothing$  100 mm).
- Position the pressure fan horizontally on the base as shown in Fig. 2.
- Optional: Fix the Prandtl pressure probe using the stand base, stand rod and Leybold multiclamp. Do not overtighten the screw of the Leybold multiclamp!
- Plug pressure sensor S,  $\pm 70$  hPa into Input A of the Sensor-CASSY 2 (Fig. 2) or into the Mobile-CASSY 2 (Fig. 4).
- Connect the 3 mm hose coming out of pressure connection  $p_1$  (top) of pressure sensor S to the 5 mm hose delivered with the Prandtl pressure probe.
- Connect the other end of the 5 mm hose to the Prandtl pressure probe outlet for  $p_{tot}$ .
- In the same way, connect pressure connection  $p_2$  (bottom) of the pressure sensor S to the  $p_s$  outlet of the Prandtl pressure probe (see Fig. 1).

*Remark: For further hints refer to instruction sheets 373 041, 373 13 and 524 066.*


**Safety notes**

Mind the safety notes in the instruction sheet of the suction and pressure fan.

Before removing the protective grid or the nozzle:

- Pull out the mains plug and
- Wait for at least 30 seconds until the suction and pressure fan comes to a complete stop.

Fig. 2: Experimental setup with Sensor-CASSY 2.




## Carrying out the experiment

*Remark: Repeat one measurement several times for estimating measuring errors.*

### a) Measuring with Sensor-CASSY 2

- If not yet installed, install the software CASSY Lab 2 and open the software.
- [Load the settings in CASSY Lab 2.](#)
- The connected pressure sensor S should be displayed if Sensor-CASSY 2 is connected to the computer via the USB port.
- Activate the connected pressure sensor S at Input A by clicking the pressure sensor S.



*Remark: Further details about connecting sensors to Sensor-CASSY 2 can be found in the CASSY Lab 2 manual or in the web help.*

- Disconnect both hoses from the Prandtl pressure probe.
- Reset the pressure sensor S by clicking  $\rightarrow 0 \leftarrow$  in the “Settings” pane when “relative pressure  $p_{A1}$ ” is marked.
- Connect both hoses to the Prandtl pressure probe again.
- Set the suction and pressure fan to its minimum speed (i.e. left limit position of fan control) and only then switch it on.
- Slowly increase the speed of the suction and pressure fan until the “relative pressure  $p_{A1}$ ” reaches approx. 1.5 hPa in a central position ( $h = 0$  cm) in front of the nozzle (distance  $x = 10$  cm.)
- Read off the “relative pressure  $p_1$ ” (in this experiment =  $\Delta p = p_d = p_{\text{tot}} - p_s$ ).
- To record the pressure values with CASSY Lab 2 click  or press the function key F9 when table “ $p_d(h)$  [autom.]” is displayed.
- Repeat these steps and read off the dynamic pressure  $p_d (= \Delta p)$  at various heights  $h$  and distances  $x$  in front of the nozzle. The grid of the nozzle may serve as a guide for the vertical position  $h$ .

*Remark: For the reset button  $\rightarrow 0 \leftarrow$  to appear in the  “Settings” pane “relative pressure  $p_{A1}$ ” has to be marked in the submenu of “CASSYs”. It is recommended to press the reset button  $\rightarrow 0 \leftarrow$  before each series of measurements.*

### b) Measuring with Mobile-CASSY

*Remark: To record the pressure values automatically, follow the instructions as described in a).*

- Switch on the Mobile-CASSY with the  key.
- Open the main menu by pressing the  key again.
- Select the submenu “Quantities” by using the arrow keys,  or , and enter the submenu using the right  key.
- Go to the submenu “p” using the right  key.
- Choose “Compensate Offset” and set the pressure to zero by pushing the right  key.
- Press  key and then the left  key to display the current pressure value.
- Set the suction and pressure fan to its minimum speed (i.e. left limit position of fan control) and only then switch it on.
- Slowly increase the speed of the suction and pressure fan until the pressure difference  $\Delta p (= p_d)$  reaches approx. 1.5 hPa in a central position ( $h = 0$  cm) in front of the nozzle (distance  $x = 10$  cm.)
- [Load the settings in CASSY Lab 2](#) and type in the pressure values in table “ $p_d(h)$  [manu.]”.
- Repeat these steps and read off the dynamic pressure  $p_d (= \Delta p)$  at various heights  $h$  and distances  $x$  in front of the nozzle. The grid of the nozzle may serve as a guide for the vertical position  $h$ .

*Remark: For further hints on using Mobile-CASSY (524 009A) refer to the instruction sheet.*


*Remark: To record more than the prepared measurement series open “Measurement” in the menu bar and select  “Append new Measurement Series”. Select table “ $p_d(h)$  [manu.]” and click  once. Open the  “Settings” pane and mark “ $p_d(h)$ ” in the submenu “Displays”. Push the button “Add new Curve” and select “ $p_d\#3$ ” in the drop down menu for “y-axis”.*

**Measuring example**

Tab. 1: Pressure profile: Dynamic pressure at a distance  $x$  from the nozzle.  $h$  corresponds to the vertical offset.

$\frac{x}{\text{cm}}$	$\frac{h}{\text{mm}}$	-100	-75	-50	-25	0	25	50	75	100
10	$\frac{p_d}{\text{Pa}}$	0	0	99	146	151	151	111	0	0
40	$\frac{p_d}{\text{Pa}}$	1	10	88	132	132	130	80	9	1

Fig. 3: Pressure profiles at two different distances:  
 $x_{\#1} = 10 \text{ cm}$ ,  $x_{\#2} = 40 \text{ cm}$ .



**Results and evaluation**

The wind speed  $v$  can be obtained with the measuring results of Tab. 1, equation (II) and the density  $\rho$  of the flow medium air:

$$\rho = 1.2 \frac{\text{kg}}{\text{m}^3}$$

From this it follows:

Tab. 2: Wind speed  $v$  calculated with the pressure values of Tab. 1.  
 $h$ : vertical offset,  $x$ : distance from the nozzle.

$\frac{x}{\text{cm}}$	$\frac{h}{\text{mm}}$	-100	-75	-50	-25	0	25	50	75	100
10	$\frac{v}{\frac{\text{m}}{\text{s}}}$	0.0	0.0	12.8	15.6	15.9	15.9	13.6	0.0	0.0
40	$\frac{v}{\frac{\text{m}}{\text{s}}}$	1.3	4.1	12.1	14.8	14.8	14.7	11.5	3.9	1.3

Within the diameter of the nozzle the pressure difference – and therefore the wind speed– varies little. Quadrupling the distance  $x$  from 10 cm to 40 cm results in a small decrease of wind speed due to the friction in the air.

In contrast to that, the lateral ( $h = \pm 75 \text{ cm}$ ) measured wind speeds are, at first, increasing with greater distance ( $x = 40 \text{ cm}$ ) from the nozzle. The cone of wind slightly fans out.

Fig. 4: Experimental setup with Mobile-CASSY.

