

## Verifying the Bernoulli equation – Measuring with the precision ma- nometer

### Objects of the experiment

- To verify that the dynamic pressure increases in reduced cross-sectional areas.
- To verify that the flow velocity increases in reduced cross-sectional areas.
- To verify that the volume flow and total pressure remain constant.





Fig. 1: Prandtl pressure probe for measuring the static pressure  $p_s$  and the total pressure  $p_{tot}$ .

Fig. 2 Bernoulli's principle: cross-sectional areas  $A_0$  and  $A_1$ , flow velocities  $v_0$  and  $v_1$ , total pressures  $p_{tot0}$  and  $p_{tot1}$ . The font size indicates the absolute value of the physical quantity.



### Principles

Bernoulli's law states the relationship between static pressure  $p_s$  and flow velocity  $v$ . The following equation applies to a friction-free, horizontally flowing stream through a stationary flow tube between two points labelled with indices 0 and 1:

$$p_{s0} + \frac{\rho}{2}v_0^2 = p_{s1} + \frac{\rho}{2}v_1^2 \quad (I)$$

Density of the air:  $\rho = 1.2 \frac{\text{kg}}{\text{m}^3}$

In particular, equation (I) states that the total pressure  $p_{tot}$  has the same value everywhere:

$$p_{tot} = p_s + p_d = \text{const.} \quad (II)$$

$p_d$ : dynamic pressure

In the experiment described here air flows through a narrowing wind tunnel. Its cross-sectional area is reducing from  $0.020 \text{ m}^2$  to  $0.015 \text{ m}^2$ . We will measure the total pressure  $p_{tot}$  and the static pressure  $p_s$  at several positions in the wind tunnel.

The flow velocities  $v_0$  and  $v_1$  at two different locations in the wind tunnel with cross-sectional areas  $A_0$  and  $A_1$  are given by the continuity equation:

$$v_0 \cdot A_0 = v_1 \cdot A_1 \quad (III)$$

The continuity equation states that the volume flow  $J = vA$  in a tube is constant. The incompressibility of air can be assumed for the occurring flow velocities in this experiment.

The definition of the dynamic pressure  $p_d = \frac{\rho}{2}v^2$  can be derived from equations (I) and (II). It allows to eliminate  $v_0$  in equation (III). By rearranging we obtain:

$$v_1 = \sqrt{\frac{2}{\rho} p_{d0}} \cdot \frac{A_0}{A_1} \quad (IV)$$

with

$$p_d = p_{tot} - p_s \quad (V)$$

The dynamic pressure  $p_d$  is determined by measuring the pressure difference with the Prandtl pressure probe. The cross-sectional area sizes are written on the ramp of the wind tunnel.

**Apparatus**

1 Wind tunnel.....	373 12
1 Prandtl pressure probe .....	373 13
1 Suction and pressure fan.....	373 041
1 Measurement trolley for wind tunnel.....	373 075
1 Precision manometer.....	373 10

*Optional:*

1 CASSY Lab 2 .....	524 220
---------------------	---------

**Safety notes**

Mind the safety notes in the instruction sheets of the wind tunnel, the suction and pressure fan and the precision manometer.

Before removing the protective grid or the nozzle:

- Pull out the mains plug and
- wait for at least 30 seconds until the suction and pressure fan comes to a complete stop.


**Setup**

Assemble the wind tunnel and the fan as shown in Fig 3. Plug the suction side of the fan into the outlet nozzle so that the air is drawn through the wind tunnel during the experiment. Ensure a clearance of approx. 1 m in front of the suction nozzle and behind the fan so that the air can be drawn into the wind tunnel without any turbulence.

- Insert the lamination filter (protection grid) into the slot at the entry gate.
- With the aid of the four screws fix the Bernoulli ramp underneath the Plexiglas hood so that the ramp height increases in the direction of flow.
- Mount the scaled sealing strip (included in the equipment for the wind tunnel) on the guiding rails.
- Carefully push the Prandtl pressure probe through the foam rubber seal of the sealing strip. Plug the Prandtl pressure probe onto the measurement trolley.
- Push the Prandtl pressure probe down so that it reaches a position approx. 2 cm higher than the highest point of the ramp.
- Slide the front of the measurement trolley to Position A of the sealing strip.
- Mount the precision manometer on the wind tunnel with the help of the stand rod, Ø 12 mm, 75 cm (delivered with the wind tunnel) as in Fig. 3.
- Connect the 5 mm hose delivered with the Prandtl pressure probe to the precision manometer's tube attachment nipple for high-pressure (left).
- Connect the other end of the 5 mm hose to the Prandtl pressure probe outlet for  $p_{tot}$  (see Fig. 1).
- In the same way, connect the tube attachment nipple for low-pressure (right) of the precision manometer to the  $p_s$  outlet of the Prandtl pressure probe (see Fig. 1).

*Remark: For further hints refer to instruction sheets 373 10, 373 12, 373 13 and 373 041.*

Fig. 3: Experimental setup with the precision manometer.



## Carrying out the experiment

### a) Measuring without CASSY Lab 2

- Set the suction and pressure fan to its minimum speed (i.e. left limit position of fan control) and only then switch it on.
- Slowly increase speed of the suction and pressure fan until the pressure difference  $\Delta p$  ( $= p_d$ ) reaches approx. 8 Pa at position A.
- Read off the dynamic pressure  $p_d$
- Note the dynamic pressure  $p_d$  together with the corresponding cross-sectional area (e.g.: 0.020 m<sup>2</sup>)
- Place the measurement trolley with the Prandtl pressure probe one position further.
- Repeat this measurement for the positions “B” to “F”.
- Repeat these steps for two more fan speeds.

### b) Measuring with CASSY Lab 2

- If not yet installed, install the software CASSY Lab 2 and open it.
- [Load the settings in CASSY Lab 2.](#)
- Slowly increase the speed of the suction and pressure fan until the pressure difference  $\Delta p$  ( $= p_d$ ) reaches approx. 8 Pa at position A.
- Read off the dynamic pressure  $p_d$ .
- Type in the pressure values and the corresponding cross-sectional area (e.g.: 0.020 m<sup>2</sup>) in table “ $p_d(A)$ ” (left side of CASSY Lab 2 window).
- Place the measurement trolley with the Prandtl pressure probe one position further.
- Repeat this measurement for the positions “B” to “F”.
- After six single measurements start a new series of measurements by clicking **#1** . Increase the fan speed a bit and put the measurement trolley back to position A.
- Repeat these steps for two more fan speeds.

*Remark: To record more than the three prepared series open “Measurement” in the menu bar and select  “Append new Measurement Series”. Then click  when table “ $p_d(A)$ ” is displayed. Open the  “Settings” pane and mark “ $p_d(A)$ ” in the submenu “Displays”. Push the button “Add new Curve” and select “ $p_d\#4$ ” in the drop down menu for “y-axis”.*

Measuring example




Fig. 4: Dynamic pressure  $p_d$  as function of the cross-sectional area  $A$ . The solid curve corresponds to a fit of a  $1/A^2$  hyperbola.

Tab. 1: Dynamic pressure  $p_d$  at positions A to F for a constant fan speed.

Position	$\frac{A}{m^2}$	$\frac{p_d}{Pa}$
A	0.020	10.4
B	0.019	12.0
C	0.018	12.8
D	0.017	14.4
E	0.016	16.1
F	0.015	19.0

Results




Fig. 5: Volume flow  $J$  as a function of the cross-sectional area  $A$  for three different wind speeds of the fan. The solid lines indicate the average value of one series of measurement.

Tab. 2: Flow velocity  $v$  and volume flow  $J$  calculated with the values of Tab. 1.

Position	$\frac{v}{m/s}$	$\frac{J}{m^3/s}$
A	4.2	0.083
B	4.5	0.085
C	4.6	0.083
D	4.9	0.083
E	5.2	0.083
F	5.6	0.084




Fig. 6: Flow velocity  $v$  as a function of cross-sectional area  $A$  for three different wind speeds of the fan. The solid curves correspond to a fit of  $1/A$  hyperbolas.

## Evaluation

The relation

$$\rho_d \propto \frac{1}{A^2}$$

derived from Bernoulli's principle and the continuity equation seems to be confirmed by Tab. 1 and Fig. 4. The best fit function is a quadratic hyperbola.

Equation (IV) seems to be confirmed by Tab. 2 and Fig. 6:

$$v \propto \frac{1}{A}$$

(Density of air:  $\rho = 1.2 \frac{\text{kg}}{\text{m}^3}$ )

The flow velocity increases with decreasing cross-sectional area  $A$ . The best fit function is a hyperbola.

The volume flow

$$J = v A$$

is constant for all cross-sectional areas (Fig. 5). The minor changes can be explained through measuring tolerance and leakages. Thus the predictions based on Bernoulli's principle and the continuity equation are verified qualitatively.

## Supplementary information

Additionally, the constant total pressure  $p_{\text{tot}}$  along the decreasing cross-sectional area  $A$  of the wind tunnel can be verified directly. Therefore only one hose is connected to the total pressure outlet (Fig. 1) of the Prandtl pressure probe as shown in Fig. 7. The hose leads to the tube attachment nipple for high-pressure (left) of the precision manometer.




Fig. 7: Experimental setup for measuring the total pressure  $p_{\text{tot}}$ .

The decreasing static pressure  $p_s$  with reducing cross-sectional areas can be verified, too. Therefore only one hose is connected to the static pressure outlet (Fig. 1) of the Prandtl pressure probe (not shown in Fig. 7). The hose leads to the tube attachment nipple for low-pressure (right) of the precision manometer.

*Remark: The decreasing static pressure  $p_s$  in reduced cross-sectional areas has already been verified in P1.8.5.4.*