
Sphärometer

Spherometer

Sphéromètre

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The spherometer serves to measure the thickness of plates and the radii of curvature of spherical surfaces, lenses etc; measuring range: 25 mm.

1. Description

The instrument consists of a stand whose three feet have steel points and are arranged in the form of a equilateral triangle. A micrometer screw with a steel point is passing through the stand in the center of the triangle. The micrometer screw carries a dial graduation of 250 divisions, ranging from 0 to 500. A vertical rule with divisions in mm is fixed on the stand.

The micrometer screw spindle has a pitch of 0.5 mm that can be read on the vertical rule. The fractions are read on the circular scale of the disk. Each scale mark corresponds to a change of altitude of 0.002 mm.

The distance between the steel points of the three feet is $d = 50$ mm, and the distance from the point of the micrometer screw to the three feet is

$$a = \frac{50 \text{ mm}}{\sqrt{3}} = 28.9 \text{ mm each.}$$

Notes

- The five-digit numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
- The specifications and illustrations are not binding in every detail for the design of the apparatus. It is our policy always to keep our manufacturing programme right up to date so that it makes full allowance for the developments acquired in all scientific and technical fields.

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$$R = \frac{e}{2} + \frac{1}{6} \frac{d^2}{e}.$$

A control of the distance is required for precise measurements.

2. Use

For measuring, place the spherometer on a plane glass plate and observe first of all the adjustment of the zero point. For this, turn down the micrometer screw so far that its pointed end is just touching the plate. If this is so, the instrument will be taken along with the rotation or it will just start to shake slightly. The friction head provided at the micrometer screw is not adapted for this precise adjustment. It serves only as protection against damages to the micrometer screw.

After having checked the zero point the measurement of the plates may be carried out in the manner described above.

For measuring the radii of curvature of spherical surfaces, place the spherometer on the spherical surface to be measured in such a way that all four steel points are uniformly touching the spherical surface. If now the reading on the micrometer screw is e and the distance between the three steel points of the equilateral triangle is d , it follows for the radius of curvature R :

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