
4/95-Kem-

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The trolley, which can move along the precision metal rail (460 81) or any flat surface, allows experiments on uniform motion and is particularly suitable as an introduction to the term velocity. The speed in both the forward and reverse direction is continuously variable.

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1 Safety note

Switch off the motor after use in order to save the batteries.
Observe polarity when inserting batteries.

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2 Description, technical data

Drive with electric motor and worm wheel

- ① Three-stage switch: forward - off - reverse
- ② Potentiometer for speed regulation:
approx. 20-60 cm/s
- ③ Battery compartment for 2 mignon cells
(685 44, not included in scope of supply)
series connected:
3 V for driving the electric motor (when
accumulators are used, the lower voltage results
in a lower top speed)
- ④ Clamp for securing metallized recording paper (337 39)

Gauge: 3.5 cm (fitting precision metal rail)
Dimensions: 12 cm 6 cm 5 cm
Mass: 175 g

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3 Sample experiments

3.1 Direct measurement of distance and time

- 3.1.1 Using a stopwatch, measure the times t required by the trolley to travel certain distances s : Result $v = s/t$
- 3.1.2 Mark the distances s the trolley travels in equal time intervals t :
Result: $s/t = v$ (const.)
- 3.1.3 Mark the distances s the trolley travels in different times t :
Result: s/t or $s/t = \text{const.}$ (= v)

3.2 Measurement of distance and time using the time recorder (337 18)

Secure metallized paper to the clip ④ of the trolley and mount on time recorder (follow the instructions). Set time marker frequency to 10 Hz. Measure the (equidistant) markings and subsequently determine the velocity.