
4/95-Di/Sf-

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Fig. 1

In time and frequency measurements on moving objects, the forked light barrier controls - either singly or in pairs - electrical timing and counting devices which are started or stopped by shorts or interruptions at their inputs.

Experiment examples:

- Track experiments
- Free-fall
- Ballistic experiments (e.g. with the large projection apparatus 336 56)
- Duration of pendulum oscillations
- Oscillation duration and frequency of leaf springs and strings

1 Safety note

Do not exceed the maximum permissible voltage (25 V DC or 15 V AC).

2 Description, technical data (cf. Fig. 1)

- ① Position mark for infrared beam
- ② Emission opening of infrared beam
- ③ Light diaphragm of receiver circuit
- ④ Light emitting diode (LED): goes out when beam interrupted
- ⑤ DIN socket, 4 pins assigned:
 - Pins 2,4: input for supply voltage
 - Pin 1: +-output for receiver circuit
 - Pin 5: --output for receiver circuit
- ⑥ Clamp for attaching apparatus to stand rods up to 14 mm dia.
- ⑦ Holding magnet
- ⑧ Stop for right-angled, tip-free mounting on rectangular profiles (Fig. 2)

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Technical data:

Required supply voltage (cf. 3.1.1):	9 V DC to 25 V DC or 6 V AC to 15 V AC
Current consumption:	110 mA
Output:	floating, electrically isolated from supply voltage

Wavelength of IR radiation:	950 nm
-----------------------------	--------

Required minimum width <i>d</i> of interrupter flag or light gap:	0.5 mm
---	--------

High-sensitivity resolution:	0.1 mm *)
Switching frequency:	max. 5 kHz

Clear fork width:	11 cm
Effective fork depth:	16 cm
Weight:	0.3 kg

*) This means:
The light barrier registers an interrupter flag or a light gap
having the width *d* with an accuracy of up to 0.1 mm.

3 Operation

3.1 Connecting timing and frequency measurement devices

3.1.1 Recommended:

*Electronic timers and frequency meters with built-in voltage
supply (6 V AC, 0.4 A) for the light barrier, e.g.*

Electronic stop-clock P	313 031
LH-Digital counter	575 40

For computer-assisted experimenting, e.g.

CASSYpack E	524 007
Software "Measuring and Evaluating"	524 112

or

Software "Motion (spoked wheel)"	524 332
in conjunction with the spoked wheel	337 461

*For correct-polarity connection of the light barrier to other appa-
ratus:*

Multi-core cable, 6-pole, 1.5 m	501 16
---------------------------------	--------

*Electronic timers and frequency meters without built-in voltage
supply for light barrier, e.g.*

Counter P	575 45
Counter S	575 47
with 12 V plug-in power supply, e.g.	562 792

Additionally required:

Voltage supply (see section 2) e.g.	
Transformer, 6 V, 12 V/30 W	562 73
Adapter cable, 4-pole, 1.5 m long	501 18

(black plug connected to voltage source, "red" to "+" and "gray"
to "⊥" of counter).

3.2 Experiment notes

Required minimum width *d* of interrupter flag or light gap:
0.5 mm (provided the accuracy is not impaired by extremely
bright room lighting or direct illumination of the receiver di-
aphragm).

In parallel or series connection of multiple forked light barriers
(337 46) or of forked light barriers and older model light barriers
(313 10), make sure all connections are made with the correct
polarity.

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