Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

7/1987 — Di/Sf —

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When used in conjunction with the aerodynamic test accessories 1 and 2 (373 07 and 08), the wind tunnel (373 12), which is connected to the suction side of the suction and pressure fan (373 04), makes it possible to carry out aerodynamic experiments in a closed working section.

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Example Experiments:

- Experimental verification of the Bernoulli equation
- Recording of the polar coordinates of an airfoil
- Aerodynamic behaviour of the air around an airfoil
- Measurement of air resistance of different bodies

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1 Safety Instructions

- Only change the base of the wind tunnel when the fan is at rest!
- Never clean the Plexiglas hood with aggressive solvents!
- Protect the guiderails ⑥ of the Plexiglas hood from dirt and scratches!
- Never use force to tighten the screws on the Plexiglas hood!
- Only move the wind tunnel in its assembled form when the floor is securely fastened and the fan is removed!

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2 Schedule of Parts Supplied, Description, Technical Data

- ① Suction nozzle, securely connected to entry gate ②
- ② Entry gate
- ③ Holes for attaching the sector dynamometer (373 14)
- ④ Recesses for the securing screws on the Plexiglas hood
- ⑤ Plexiglas hood
- ⑥ Two polished guiderails for mounting the dynamometer car (from aerodynamic test accessories 1, 373 07)
- ⑦ Slot for the airfoil holding bars and the pressure head
- ⑧ Four knurled screws for securing the Plexiglas hood ⑤ to the entry and exit gates ② and ⑨
- ⑨ Exit gate with recesses, same as ④, with fitted safety filter
- ⑩ Screw socket for screwing in the stand rod ⑳
- ⑪ Diffuser, secured to the exit gate ⑨, for attachment of the suction and pressure fan (373 04)
- ⑫ Level floor with expanded foam seal all around for airtight sealing of the wind tunnel floor during aerodynamic tests
- ⑬ Four securing screws for screwing on the floors ⑫ or ⑭ on the entry and exit gates
- The floor can be exchanged for Bernoulli ramp ⑮ or Venturi floor (373 21)
- ⑭ Slot for inserting the filter ⑮
- ⑮ Filter for inserting into the slot ⑭ of the entry gate for protecting the working section from foreign matter drawn in by suction and for lamination of the air current
- ⑯ Bernoulli ramp, used in place of the level floor ⑫ during experiments on Bernoulli's Law;
Scale: cross-sections of the experimental area in the wind tunnel
- ⑰ Sealing strip for sealing the slot ⑦ in the Plexiglas hood during experiments with the Bernoulli ramp ⑯, with position marks for the pressure head (373 13) placed on the dynamometer car. The marks correspond to those on the cross-section scale on the Bernoulli ramp ⑯
- ⑱ Transparent plastic hood for sealing slot ⑦ in the Plexiglas hood during experiments on air resistance


Fig. 1

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- ⑲ Back wall made of black plastic with horizontal white orienting lines. With two holes for mounting behind the Plexiglas hood (see Section 3.1)
- ⑳ Stand rod, diameter 12 mm, with screw thread M8, length 75 cm, for screw socket ⑩ ; for securing the precision manometer (373 10)

Not shown:

Dust cover, fits over the assembled wind tunnel, the suction nozzle and the diffuser

Dimensions (without fan): 36 cm x 42 cm x 113 cm
Weight: 6 kg

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3 Operation

3.1 Assembly (Fig. 2)

Place suction nozzle ① on a table with the suction opening facing downwards.

Unscrew all knurled screws ⑧ of the Plexiglas hood ⑤ so that the thread no longer projects inside the hood.

Orientate Plexiglas hood so that the asymmetrical slit ⑦ is positioned exactly on top as shown in Fig. 2; then push hood vertically into the entry gate ② as far as possible (Fig. 2) so that the threads of the knurled screws ⑧ pass into the recess ④. Tighten both knurled screws.

Place exit gate ⑨ over the upper end of the Plexiglas hood in a similar way and secure by tightening the screws ⑧.

Close the remaining opening of the Plexiglas hood with the level floor ⑫ or the Bernoulli ramp ⑬. To do this, screw the four screws ⑩ in the corners of the appropriate floor into the corresponding screw-sockets of the entry and exit gates.

Place the wind tunnel on a table with the suction nozzle ① facing left.

If the back wall ⑲ is to be attached: unscrew the two rear knurled screws ⑧ fully, push through the holes in the back wall and retighten so that the back wall is fixed behind the Plexiglas hood.

Place the fan on its base and, taking note of the direction of flow, push into the diffuser in such a way that the air is sucked through the wind tunnel during the experiment.


Fig. 2

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Fig. 3 Set-up for experiments on Bernoulli's Law: Measurement of dynamic pressure as a function of cross-section; wind tunnel fitted with Bernoulli ramp and sealing strip.

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Fig. 4 Set-up for aerodynamic experiments in the wind tunnel with level floor: Measurement of lift and air resistance on an airfoil with different angles of incidence (recording of the polar coordinates).


Fig. 4

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4 Maintenance, Storage

Store the wind tunnel under the dust cover supplied in order to prevent the Plexiglas hood from gathering dust and yellowing.

Only clean with non-alcoholic, non-aggressive rinsing agents and water or special cleaners (e. g. Plexiklar®)!