
6/95-Sf-

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## Lampengehäuse Lamp Housing

Fig. 1

The lamp housing (450 60) with incandescent lamp 6 V, 30 W (450 51) is used as a source of light in many optical experiments as well as for projections, production of luminous pointers on mirror instruments and similar. The advantages of this source of light are high luminous intensity of the filament, easy handling and simple adjustability.

## 1 Safety Notes

Max. permissible voltage for lamp 6 V, 30 W:  
8 V (for short periods)

Do not cover the venting holes at the rear of the lamp housing.

The lamp housing will become very hot!

## 2 Description

- ① Lamp insert, turnable and shiftable in lamp housing ⑤
- ② 3 adjusting screws for lamp socket ④
- ③ Venting holes
- ④ Lamp socket E 14 for lamp 6 V, 30 W (450 51)
- ⑤ Lamp housing
- ⑥ Holder for condenser (460 20 or 17)
- ⑦ Stand rod, 11 cm x 1 cm dia.
- ⑧ Connecting lead with two 4-mm plugs

Dimensions of housing (unextended):  
approx. 12 cm x 7 cm dia.

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Transformer 2 V to 12 V, 10 A 562 75

Condenser, e.g.

Aspherical condenser (with  
 diaphragm holder, diaphragms,  
 imaging objects, luminous-  
 pointer slide) 460 20

or

Single-lens condenser 460 17  
 (with sliding aperture selector)

Stand material, e.g.

Saddle base 300 11  
 or (for arrangements with  
 optical axis)  
 Optical bench, e.g. 460 43  
 Leybold multiclamp 301 01  
 Large stand base 300 01

For screwing in the incandescent lamp, remove the lamp insert  
 ① from the lamp housing ⑤. Align the lamp by means of ad-  
 justing screws ② so that its axis is perpendicular to the base of  
 the lamp insert.

Depending on experimental requirements, turn lamp inset ① to  
 adjust either to the horizontal direction (e.g. for experiments  
 with grazing incidence of light on a screen, according to  
 Fig. 2.1) or perpendicular direction (e.g. for illuminating a vertical  
 slit, see Fig. 2.2).

For experiments on the optical bench, fit stand rod ⑦ in the  
 multiclamp so that it is flush with the lower edge of the clamp  
 (coarse adjustment of axis in connection with optical elements  
 460 01 and following numbers, which are similarly fastened in  
 the multiclamp, see Fig. 2.2).

For focusing the light, insert condenser in holder ⑥ and shift  
 lamp insert ① in lamp housing ⑤ until the desired adjustment is  
 obtained (e.g. parallel light, see Fig. 2.1).


Fig. 2.1

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