
5/95-Di/Sf-

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Asphärischer Kondensator Aspherical Condenser

Fig. 1

The condenser is used for focusing the light of the incandescent lamp 6 V, 30 W (450 51) in lamp housing (450 60). The plug-on diaphragm holder can hold framed slides, diaphragms, image objects and filters having the dimensions of a 5 cm x 5 cm slide frame (the scope of supply includes 5 diaphragms and image objects).

1 Description, scope of supply, technical data

- ① Condenser in socket, 4 cm x 7.5 cm dia.
 - (1.1) Aspherical condenser, 60 mm dia., $f = 49$ mm
 - (1.2) 4-mm holes for attaching the diaphragm holder ②
 - (1.3) Muff for fitting into the lamp housing (450 60)
- ② Diaphragm holder, light aperture 4.5 cm x 4.5 cm
 - (2.1) Clamping rail for mounting diaphragms etc. up to a thickness of 3 mm
 - (2.2) 4-mm pins, fitting into the holes (1.2)
- ③ to ⑦ Diaphragms and image objects, 5 cm x 5 cm
- ③, ④ Diaphragms with holes,
opening: 12 mm dia.
6 mm dia.
- ⑤ Slit diaphragm,
slit width: 1 mm
- ⑥ Arrow diaphragm (image object);
dimensions of arrow: 1 cm x 0.2 cm
- ⑦ 'Luminous pointer' slide with division mark (width 0.5 mm) on transparent field (dia. 20 mm), for producing a 'luminous pointer' in mirror-scale arrangements


Fig. 2.1


Fig. 2.2

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2 Operation

2.1 Focusing the light beam of the incandescent lamp 6 V, 30 W (450 51) in lamp housing (450 60)

To produce divergent, parallel and convergent light beams, move the lamp insert ① in longitudinal direction. Beam control is either by means of a screen according to Fig. 2.1 or, with horizontal lamp filament, by means of a screen or sheet of paper according to Fig. 2.2. If necessary, adjust the lamp in the optical axis by turning the screws ②.

2.2 Luminous pointer (e.g. for experiments with the centrifugal force apparatus, 347 22)

To measure the deflection of a mirror in a luminous-pointer arrangement, mount lamp housing with incandescent lamp (450 60/51) into the experimental assembly using stand equipment, so that the light falls onto the mirror.

Push the 'luminous pointer' slide ⑦ into the holder ② so that the division mark has the same direction as the graduation on the scale onto which image is projected. Adjust lamp filament accordingly.

Produce an image of the 'luminous pointer' slide on the mirror by shifting the lamp insert, and then move the complete lamp housing until the projected mark is sharply focused onto the scale.