
6/88 — lbg —

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This model is used to explain and demonstrate by experiment the guidance of light by total reflection. This effect is industrially used in fiber optics.

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1. Description

The apparatus consists of a helically shaped rod of acrylic glass, 500 mm in length and 10 mm in diameter, with two windings of approx. 50 mm dia. The acrylic glass rod is fitted on a rider for mounting to an optical bench. The light incident on one face of the rod is guided through the rod and its windings by repeated total reflection and is emitted at the other face of the rod as a divergent light beam. Due to the unavoidable inhomogeneities of the acrylic glass and to the somewhat small winding diameter as compared with the diameter of the rod, the total reflection is slightly disturbed. A small fraction of the light guided through the interior of the rod is emitted tangentially in the direction of the light, visibly for the observer, along the bent outer surface of the helix.

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2. Carrying out the Experiment

For an off-hand experiment the light guide model can e. g. be placed on a stand base and one face of the acrylic glass rod be illuminated by a high-duty pocket lamp with convex lens. The light guided through the rod and emitted at the other end is then demonstrated by a translucent screen, e. g. a piece of tissue paper.

The efficiency of the experiment can be improved by mounting the assembly carefully on an optical bench using a high-luminosity focusing lens system.

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
For the experiment, lamp house with lamp are fitted to one of the riders of the optical bench. With inserted condenser and by moving the lamp house, the filament is formed on an approx. 20 cm distant focal plane; thereafter the condenser is provided with the pin diaphragm. The light guide model is then arranged on the second rider in front of the lamp with the condenser in such a way that the light emitted from the diaphragm aperture is centrally incident on the face of the helical plastic rod. This entrance area should then be about 10 mm before the diaphragm aperture. If necessary, the centering of light on the face of the plastic rod must be improved by adjusting the lamp house, possibly also by readjusting its set-screws. At last, the pull-out of the lamp house is adjusted until the light emitted on the other face of the light-guide model at an angle of 90° to the direction of the optical bench shows optimum brightness.

The light emitted from the sides of the light-guide model is experimentally demonstrated by keeping the translucent screen off-hand into the deflected light beam in the intermediate vicinity of the light-guide model. Due to the divergence of the emitted light the brightness on the illuminated translucent screen is reducing as the distance from the light-guide model increases.

When facing the light emitted from the light-guide model on the optical bench, the lower sides of the windings in the helical acrylic glass rod seem to be distinctly brighter than the other areas. Viewing however in direction of the light parallel to the optical bench over the lamp house, the same phenomenon can be observed on the upper sides of the helical windings. A minor part of the light guided in the acrylic glass rod is emitted from the helical windings of the rod tangentially to the direction of the light guided in its interior. The total reflection is slightly disturbed as with the internal reflection along the helical rod it cannot be strictly maintained and, in addition, the acrylic glass shows traces of inhomogeneities which cannot be avoided in fabrication.

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