01/89 -Kp-

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Fresnel's mirrors can be used to observe interference fringes which confirm the wave theory of light. It is possible to determine quantitatively the wavelength of the light source and the angle of incidence between the two planes of the Fresnel's mirrors. This experiment can be set up and adjusted quickly and easily, either for demonstration purposes or for use in student labs.

The interference fringes are created when the light generated by source A (can be regarded as point source) is reflected from the two mirrors, the one being raised slightly from the plane to form a shallow angle. Coherent waves are propagated from the two virtual sources of light thus simulated (Fig. 1).

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Fig. 1
Schematic representation of the light propagation in the experiment with Fresnel's mirrors

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At large distances L (see Fig. 2), the following formula can be used to calculate the spacing between two adjacent maxima (crests) or minima (valleys).

$$\frac{d}{L} = \frac{\lambda}{a} \quad \lambda = \frac{a \cdot d}{L}$$

whereby
d: distance between two adjacent maxima (or minima)
L: distance between light sources and projection plane
 λ : wavelength of the light
a: distance between the two virtual light sources

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A simple optical configuration is used to determine distance a between the two virtual light sources. To do this, a lens is used to project the two virtual light sources in sharp focus on an observation screen (Fig. 3):

$$\text{Here either } \frac{a}{g} = \frac{B}{b} \quad \text{or} \quad a = \frac{B \cdot g}{b} \text{ applies,}$$

whereby

g : distance between the lens and the light source

b : distance between the lens and the observation screen

B : distance between the two projected images of the light sources

Fig. 3

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Fig. 3

Determining distance a between the two virtual light sources A' on the basis of the equation for the projection

The following formulae are used to determine the included angle α of inclination between the two mirrors (see Bergmann/Schäfer, Lehrbuch der Experimentalphysik, Optik)

$$\alpha = \frac{a}{2 \cdot R} \quad \text{or} \quad \alpha = \frac{\lambda \cdot L}{2 \cdot R \cdot d}$$

whereby

α : angle of inclination between the two mirrors, in radians

R : distance between the light source (slit) and the Fresnel's mirrors

Where the angle of inclination is known, the wavelength can be established without having to determine distance "a" between the two virtual light sources (see above); this is done as follows:

$$\lambda = \frac{2 \cdot \alpha \cdot R \cdot d}{L}$$

1 Description and technical data

Fig. 4

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Fig. 4

Fresnel's mirrors on backing plate (471 04)

Two mirrors inclined toward each other ① and glued to a metal plate ②.

Dimensions of the individual mirrors: 30 mm x 40 mm

Angle of inclination: approx. 10°

Dimensions of the metal plate: 50 mm x 100 mm

2 Setup and adjustment

2.1 Setup with STM*- equipment (preferably for student laboratory use)

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1 Mirror box with stand rod	459 17
1 Measuring magnifier	459 19
1 Red filter	from 467 95
1 Steel tape measure, 2 m	311 77

The two rails (initially without lens H) are to be set up as shown in Fig. 5, at a slight angle one to the other, so that the Fresnel's mirrors are positioned in the beam of light, about 90 cm from the slit. Adjust the Fresnel's mirrors so that the reflected images of the slit can be seen in the mirror box without using the measuring magnifier.

To observe the interference fringes, lay the measuring magnifier in the mirror box and reduce the size of the slit to about 0.1 mm.

To determine the distance between the virtual slits, position lens H ($f=30$ cm) between the Fresnel's mirrors and the mirror box in such a way that both slits are in sharp focus.

This experiment can be conducted in normal room light.

2.2 Setup with demonstration equipment

Additional equipment required:

1 Simple bench clamp	e.g. 301 07
1 Small optical bench	460 43
5 Leybold multi-clamps	301 01
1 Halogen lamp, 12 V, 50/100 W	450 64
1 Transformer, 12 V, 100 W	e.g. 562 75
2 Connecting leads	e.g. 501 28
1 Adjustable slit	460 14
1 Prism table with holder	460 25
1 Lens with holder, $f=200$ mm	460 04
1 Holder with spring clips	460 22
1 Monochromatic filter	e.g. 468 09
1 Saddle base	e.g. 300 11
1 Translucent screen	441 53
1 Steel tape measure, 2 m	e.g. 311 77
1 Precision vernier calipers	e.g. 311 54


Set up the configuration as shown in Fig. 6 (initially without the $f=200$ mm lens). Set the adjustable slit horizontally, so that it is parallel to both the lamp filament and the mirror. Set the width of slit at about 0.1 mm. Adjust the height of the prism table so that the interference fringes are seen on the screen.

In order to determine the distance between the virtual light sources position the lens ($f=200$ mm) so that a sharp image of both light sources may be observed; the distance B between them can then be measured using the vernier calipers.

Conduct this experiment in a darkened room.

* Science Teaching Modules

Fig.5 Configuration with STM equipment

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Fig 6
Configuration with demonstration equipment

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