11/93-Kem-

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This device is for demonstrating the linear electro-optical effect (Pockels' effect), whereby the double refraction of a crystal changes under the influence of an electric field. In the case of lithium niobate, this effect is proportional to the applied voltage. The optic axis of the crystal lies parallel to the electric field of a plate capacitor and perpendicular to the principal axis of the geometrical optical path (transverse configuration). The experiments are performed with a polarised laser (471 840).

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1 Safety notes

Do not exceed a maximum voltage of 2kV.
Avoid using voltages which present shock hazards.
When working with the laser make sure that all of the safety measures provided in the appropriate instruction sheet have been met.

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2 Description, technical data

The cuboid crystal lies in a small plate capacitor. This configuration - which can be rotated around the axis of the optical path - is integrated into a holder on a rod.


Fig. 2

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Electrode configuration

- ① Electrodes of plate capacitor
- ② Crystal
- ③ Direction of field and optic axis of crystal
- ④ Optical path

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
Crystal: Lithium niobate (LiNbO₃)
 Dimensions of crystal: 2 mm x 3 mm x 20 mm
 Optical path length: 20 mm
 Plate capacitor
 Distance between plates: 2 mm
 Surface area of each plate: 3 mm x 20 mm

Half-wave voltage: approx. 380 V
 Angular range: ± 90°

Electrical connections: 4-mm safety sockets

Diameter of holder: 130 mm
 Dimensions of rod: 85 mm x 10 mmØ

3 Operation

Additionally required:
 1 He-Ne laser, linearly polarized 471 840
 1 Polarising filterfrom 472 40

1 High-voltage power supply, 10 kV 522 37
 1 High-voltage cable 501 05
 Cables

1 Optical rail, e.g.
 Optical bench with standardized V-shaped cross-section 460 32
 1 Rider 460 352
 or
 1 small optical bench 460 43
 Leybold multiclamps (used as riders) 301 01

For additional equipment, refer to sample experiments in Section 4.
 When setting up the experiments, always arrange the crystal at an angle of 45° to the polarisation plane of the laser, aligning it so that no reflections occur at the internal sides of the crystal. Position the polarising filter either parallel or perpendicular to the laser's polarisation plane.

4 Sample experiments

4.1 Demonstrating the electro-optical effect (Observation on a screen)

Additionally required:
 1 Lens in frame, $f = 5$ mm 460 01
 1 Lens in frame, $f = 50$ mm 460 02
 1 Translucent screen 441 53

Fig. 3

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Optical arrangement for demonstration of double refraction
 ① Laser
 ② Diverging lens, $f = 5$ mm
 ③ Lens, $f = 50$ mm
 ④ Pockels' cell
 ⑤ Polarising filter
 ⑥ Screen

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To demonstrate the induced double refraction, assemble the equipment as shown in Fig. 3. Position the lens ($f=50$ mm) and Pockels' cell so that light is focussed in the crystal. If necessary, rotate Pockels' cell slightly around its own axis until hyperbolic lines (due to natural double refraction) are obtained on the screen. When a voltage is applied, this lined pattern changes (induced double refraction). Each time the voltage is increased by the so-called half-wave voltage, the image is inverted, i.e. bright and dark areas reverse positions (confirming the linear nature of the effect).

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4.2 Determining the half-wave voltage (with a photoelement)

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- ① Laser
- ② Pockelszelle
- ③ Polarisationsfilter
- ④ Fotoelement.

Optical arrangement for determining the half-wave voltage

- ① Laser
- ② Pockels' cell
- ③ Polarising filter
- ④ Photoelement

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Assemble the equipment as shown in Fig. 4. Increasing the voltage causes the short-circuit current (max. approx. 100μ A) or the no-load voltage of the photoelement to change. The half-wave voltage is given by the difference between the voltage applied across the Pockels' cell at a minimum and maximum value of the received signal.

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4.3 Light modulation for transmission of sound frequencies

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AC/DC amplifier:

100-fold to 1000-fold, AC

Assemble the equipment (Fig. 4) and electrical circuit as shown in Fig. 5. Connect the photoelement to the AC/DC amplifier. Apply the amplified signal to the speaker. Vary the voltage until a clearly audible signal is obtained.