
Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

11/96-Pr/Sf-

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This device supplies a DC voltage which is continuously adjustable up to 10 kV. The output current is limited so that the voltage is non-hazardous (contact-safe). The output current is less than 100  $\mu$ A at voltages over 5 kV.

The voltage is set either by means of the built-in potentiometer or externally via an extra-low voltage control input. An integrated voltmeter displays the output voltage in digital form.

#### Examples for use:

- Electrostatics experiments
- Experiments on radioactivity (operating a Geiger counter, 546 28, and the ionization chamber, 546 25)
- Operating electron beam and gas discharge tubes, e.g.:
  - Electron beam tubes (555 10/11/12/17/214/216/218/220/234);  
Field emission microscope (554 60)
  - Spectral tubes (467 64 - 69)
  - Discharge tube P (555 222); discharge tube, open (554 16);  
Hittorf's tube (554 36)

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## 1 Safety notes

### Read this Instruction Sheet carefully!

- When this device is used according to this Instruction Sheet, experimenting with high voltages represents no safety hazard!
- This device supplies contact-safe high voltage.  
In accordance with EN 61010-1 (VDE 0411), a part is deemed to be contact-safe when, at voltages greater than extra-low voltage (60 V DC voltage),  
the current through an induction-free resistance of 2 k $\Omega$  is not greater than 2 mA for DC  
and additionally, the charge for voltages up to 15 kV is less than 45 F128MmC.
- Do not connect capacitors with a capacitance  $\geq 2.5$  nF (at 10 kV), as according to VDE 0411 a danger of hazardous contact exists at 10 kV for capacitances above 4.5 nF (approx. 2 nF already exists in the high voltage power supply).
- Do not connect multiple power supplies in series.
- Always switch off the power supply before altering the experiment setup.
- Before switching on the device, turn potentiometer ② all the way to the left (output voltage 0).
- Connect only resistors in plastic housings suitable for high voltage (e.g. 536 25) to the outputs; do not use resistors in metal jackets (old design) or mounted in STE-housings (flashover!).
- Take the ground connection of socket ⑧ and the bottom socket of input ⑨ into consideration.

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- Use high voltage cables (501 05), as the insulation of connecting leads designed for extra-low and low voltages are generally not high-voltage proof.  
If high voltage cables are not available, keep a minimum distance of 4 cm at 10 kV (or a correspondingly lesser distance for lower voltages) between the connecting leads and conducting surfaces (benchtop, experiment apparatus) in order to prevent high-voltage flashovers.  
You may need to incline the device on its folding feet to achieve a greater distance between the high voltage outputs and the benchtop.
- Set up the experiment so that it is not possible to touch either non-insulated parts or cables and plugs inadvertently.
- **Operating electron beam and gas discharge tubes:**  
As evacuated tubes can emit x-rays, these are subject in Germany to the applicable x-ray control regulation when operated at a voltage  $U \geq 5$  kV; this regulation defines such arrangements as noise radiators. In Germany, such noise radiators may not exceed the legally mandated limit values for the dose rate ( $\leq 1 \mu\text{Sv/h}$  at a distance of 0.1 m).  
As the current is internally limited to max. 100  $\mu\text{A}$  at voltages  $\geq 5$  kV, the high voltage power supply, 10 kV makes an important contribution to preventing inadmissible x-ray radiation.  
In the case of some tubes, you will need to check whether the tube manufacturer permits the operation of the tube with the electrical data specified for the high voltage power supply 10 kV (this condition is fulfilled for the tubes listed on page 1).  
X-ray protection requirements may vary from country to country. If you are unsure about the regulations which apply to you, please consult your local authorities.
- If the mains connection voltage specified on the rating plate (rear of device) does not match your local power mains, return the device suitably packed for transport to Leybold Didactic for conversion.


Fig. 1

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## 2 Description, technical data, scope of supply

- ① On-off switch with operating indicator lamp
- ② Potentiometer for direct, continuous setting of the output voltage or for preselecting a maximum value for external control of the output voltage via input ⑨
- ③ 2 1/2-digit display for output voltage  
Measuring range: 0 to 10 kV  
Accuracy of voltmeter: 3 % of end value
- ④ Changeover switch for activating outputs
  - Right position: output ⑥
  - Left position: output ⑦
  - Middle position: outputs ⑥ and ⑦ connected in series
- ⑤ LEDs for indicating the active outputs ⑥, ⑦ or ⑥ and ⑦ in series

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- ⑥, ⑦ Outputs (two pairs of 4-mm safety sockets) for tapping the voltage set using potentiometer ② or controlled externally via input ⑨;  
output voltage current-limited and floating-ground
- Output ⑥: 0 to approx. 5 kV  
Tappable current: max. 2 mA (short-circuit current)
  - Output ⑦: 0 to approx. 5 kV  
Tappable current: max. 200 µA (short-circuit current)
  - Outputs ⑥ and ⑦  
connected in series: 0 to approx. 10 kV resp. -5 kV...0...+5 kV at grounded center tap  
Tappable current: up to 5 kV max. 200 µA (short-circuit current);  
> 5 kV max. 100 µA (short-circuit current)
- Residual ripple < 1 V<sub>pp</sub>  
Outputs activated using changeover switch ④;  
active input indicated via LED (5.1) or 5.2 or (5.1) and (5.2)  
Outputs voltage-free in max. 30 s after switch-off
- ⑧ Ground socket, electrically connected to protective earth
- ⑨ Control input (4-mm safety sockets) for external control of output voltage within the upper limit set with potentiometer ② (see section 3.2), bottom socket connected to protective earth; output voltage controlled via external resistor or external extra-low voltage
- Control voltage: 0 to 5 V DC or  
0 to 5 V<sub>p</sub>, max. 1 Hz
- Corresponding  
output voltage: 0 to approx. 5 kV (output ⑥ or ⑦)  
0 to approx. 10 kV (outputs ⑥ and ⑦ in series).
- ⑩ (Rear of housing)  
6.3 AC output (4-mm safety sockets), maximum load capacity 2 A, high-voltage proof up to 10 kV for heating-voltage supply, e.g. of cathode ray tubes.

Appliance-plug connection with integrated fuse holder for primary and reserve fuse also on rear panel.  
Supplied with mains power lead.  
Two folding feet in bottom for inclining the device.

#### Additional technical data

Mains connection voltage: 230 V AC; 50/60 Hz or 115 V AC (according to rating plate on rear of device)

Power consumption: 30 VA

Primary fuse: Fuse T 0.315 for 230 V AC  
T 0.63 for 115 V AC

Dimensions: 20 cm x 21 cm x 23 cm

Weight: 2.5 kg

### 3 Operation

#### 3.1 General

Important: observe all safety notes (section 1)!  
Always make sure the power supply is switched off and potentiometer ② is in the zero position (turned all the way to the left) before setting up the experiment.  
Ground your high voltage circuit at a suitable point, and avoid ground loops.  
When connecting a gas discharge tube (e.g. spectral tubes 467 64-69 or discharge tubes 554 16, 554 36, 555 222) connect a 100-kΩ resistor (536 25) in series to prevent oscillations in the tubes due to firing phenomena.

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
Select the output according to the experiment setup, e.g.:

- Electrostatics experiments: output ⑦ (0 to approx. 5 kV, 100  $\mu$ A) or outputs ⑥ and ⑦ in series (0 to approx. 10 kV; 200  $\mu$ A to 5 kV, above this 100  $\mu$ A)
- Operating ionization chambers and gas discharge tubes: Output ⑥ (0 to approx. 5 kV, 2 mA) In exceptional cases, connect outputs ⑥ and ⑦ in series (0 to approx. 10 kV), e.g. for discharge tube (554 16) and the field emission microscope (554 60); refer to the respective Instruction Sheet.

### 3.2 Externally controlling the output voltage

Recommended e.g. for experiment setups

- which are sensitive to electrical fields caused by the experimenter ("remote control" of high voltage; see Fig. 2, Fig. 3.1/2)
- or
- which require fine adjustment of the high voltage (see Fig. 4)
- or
- in which the high voltage is computer controlled (only possible in experiments without spark gap).

Fig. 2

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Switching the output voltage on and off via an external switch ("remote control")

Connect switch (a) (e.g. 505 45) via long connecting leads; with switch (a) initially open, preselect the required maximum output voltage on potentiometer (2). Switch on the output voltage by closing switch (a), and deactivate it by opening switch (a).


Fig. 3.1

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Continuously varying the output voltage via a rheostat at control input (9) Set potentiometer (2) to the full right deflection and connect rheostat R to input (9) via long connecting leads; the value of the rheostat is determined by the maximum adjustable output voltage  $U_{max}$ :

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Continuously varying the output voltage via an externally adjustable control voltage (0...5 V) at input (9)

Either set potentiometer (2) to full right deflection or, with input (9) open, set it to the desired maximum output voltage value.


Fig. 4

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Detecting  $\alpha$  and  $\beta$  radiation using a Geiger counter (546 28) and oscilloscope; distinguishing the types of radiation by means of pulse amplitude analysis.

Fine adjustment of the high voltage in an extremely narrow proportional range of the Geiger counter using a ten-turn potentiometer (577 93).  
 With input ⑨ open, use the potentiometer ② to set an upper limit for the output voltage of approx. 3 kV.

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#### 4 Replacing the primary fuse

Pry out insert ① with holder for primary fuse ② and spare fuse ③ (Fig. 5.1).

Replace the defective fuse ② using a new fuse ③ which has been checked for the correct rating (Fig. 5.2).

Insert the new spare fuse ③ (see section 2 for rating) and put the insert ① back in the device.


Fig. 5.1


Fig. 5.2