Elektrometerverstärker
Electrometer Amplifier

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The electrometer amplifier (gain factor 1) is an impedance transformer with extremely high-resistance voltage input ($\geq 10^{13} \Omega$) and low-resistance output ($\leq 1 \Omega$).

The unit is used with a connected voltmeter (e.g. indicator instrument, recorder) to conduct the following measurements:

- 1.1 For quasi-static voltage measurement up to 10 V DC.
- 1.2 For high-resistance measurement of voltages greater than 10 V with an ohmic voltage divider consisting of STE resistors.
- 1.3 For quasi-static measurements of voltages greater than 10 V DC with a capacitive voltage divider consisting of STE capacitors.
2. For measurement of very small currents (1 V at output $\approx 10^{-6}$ A to 10^{-10} A), using STE plug-in resistors from M Ω to 10 G Ω as shunts.
3. As a charge meter (1 V at output $\approx 10^{-6}$ As to 10^{-9} As), with STE capacitors from 10³ nF to 1 nF.

Examples of experiments:

- Measuring charges and voltages in electrostatics
- Measuring charges and voltages on a plate capacitor (544 22 or 23); determination of ϵ_0 .
- Voltage measurement in liquid conductors; electrochemical investigations with a single-rod electrode; recording of equipotential lines in an electrolytic tank (545 09).
- Air ionization due to combustion gases or α radiation
- External photoelectric effect (Hallwachs effect).
- Determination of h with the photo-cell (558 77)
- Franck-Hertz experiment

Literature:

Book "Electrostatics 2: Charges and Fields"
student's experiments with information for the teacher (589 172)

1 Safety Notes

- This device is intended for use in circuits which are operated using protective extra-low voltage:
 - Maximum permissible DC voltage: 60 V
 - Maximum permissible AC voltage: 25 V
- Higher voltages may only be used when it is ensured that they are immediately switched to lower voltages in the case of contact with live (voltage-carrying) parts. This is ensured in the circuits described in section 3.1.
- Max. permissible AC supply voltage: 25 V AC; optimum: 12 V AC
- Do not apply any external voltages to the output!
- Voltage divider circuits for measuring voltages greater than 10 V (see section 3.2.2/3 and fig. 2.2/3) must be fitted with R_2 resistors or C_2 capacitors respectively, with sufficient dielectric strength for the applied voltage; comply with equipment recommendation in sections 3.1.2 and 3.1.3; in particular, be sure to use capacitors which are labelled "650 V"!

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Do not clean housing with aggressive solvents (just wipe with a damp cloth).

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- ① High-resistance amplifier input
Socket (1.1) together with sockets connected to it internally as shown in circuit diagram printed on housing are connected to earth as reference potential;
Socket (1.2) together with sockets connected internally or externally are for connection of the positive or negative potential of the input voltage (with reference to earth)
Socket spacing: 19 mm
(for plugging in STE resistors and capacitors)
- ② Auxiliary input (socket spacing 19 mm), e.g. for plugging in a voltage divider resistor or capacitor, a Faraday's cup (546 12) or an open ionization chamber (see fig. 5/6).
- ③ Low-resistance amplifier output for connection of any voltmeter (max. range 10 V) requiring less than 5 mA current e.g. indicator instruments, XY or YT recorders.
- ④ Trimmer for offset voltage (operation generally not required; see section 3.4).
- ⑤ Connecting sockets for AC supply voltage (12 V AC), from which the DC voltage required for the operational amplifier is derived by internal rectification.

Inputs (5.1) and (5.2) are connected in parallel and can be used as desired.

- (5.1) 4-mm socket pair for connecting the voltage via connecting leads with 4-mm plugs
- (5.2) Socket for cannon plug of plug-in supply unit 12 V AC/20 W (562 791 for 230 V AC or 662 792 for 115 V AC).

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Technical data

Gain:	1.00
Input resistance:	$\geq 10^{13} \Omega$
Input current:	$\leq 0.5 \text{ pA}$
Internal input capacitance	$\leq 50 \text{ pF}$
Input voltage range:	$\pm 10 \text{ V DC}$ 1 kV DC (from power supplies) 10 kV DC (generated using friction rods)
Overvoltage protection:	
Output voltage excursion:	$\pm 10 \text{ V DC}$ for a supply voltage of 12 V AC
Output current:	5 mA, proof against constant short circuit
Output resistance:	$\leq 1 \Omega$
Dimensions:	approx. 11.5 cm x 11.5 cm x 3 cm
Weight:	0.15 kg

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3 Operation

Important!

Because of the high sensitivity of the electrometer amplifier, electrostatic induction effects may cause interference in some experiments (e.g. charge measurement with a Faraday's cup; measurement of ionization currents; Franck-Hertz experiment).

Remedy: connect the experimenter via a connecting rod (532 16) and experiment cable to an earth socket of the amplifier; ground the earth connection.

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3.1 Additional requirements

Voltage source 12 V AC for supply to electrometer amplifier, e.g.

Plug-in supply unit 12 V~/20 W with cannon plug (562 791 for 230 V AC or 562 795 for 115 V AC)

or

6 V/12 V transformer; 30 W 562 73

or

Power supply 450 V DC 522 27

Multimeter for DC voltages $\leq 10 \text{ V DC}$, e.g.

AV-meter 531 94

or

Multimeter METRAMax 2 531 100

or

Recorder e.g. 575 662

DC voltage source for the experiment circuit in which voltage, current or charge is to be measured, e.g.

Low voltage source e.g. 521 35

or

for non-hazardous voltages $\geq 42 \text{ V}$:

Power supply 450 V DC 522 27

or

High-voltage power supply unit 10 kV 521 70

or pair of friction rods 541 00

Connection rod 532 16

STE high-resistance resistors and STE capacitors: see sections 3.1 to 3.3.

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3.2 Voltage measurement

3.2.1 Quasi-static measurement of voltages up to 10 V DC (see fig. 2.1)


Fig. 2.2

Hochohmige Messung von Spannungen über 10 V
High-impedance measurement of voltages greater than 10 V

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3.2.2 High-resistance measurement of voltages greater than 10 V DC

Principle (see fig. 2.2):

The total voltage to be measured

$$U_0 = U_1 + U_2 \quad (1)$$

is divided by an ohmic voltage divider R_1, R_2 ; U_1 is the
amplifier input voltage U_e . The following applies:

$$\frac{U_1}{U_2} = \frac{R_1}{R_2} \quad (2)$$

or

$$U_1 = \frac{R_1}{R_2} \cdot U_2 \quad (3.1)$$

and

$$U_2 = \frac{R_2}{R_1} \cdot U_1 \quad (3.2)$$

If R_1 is low compared to R_2 and therefore low compared to U_2
due to (2), U_1 may be neglected in (1). Consequently, from (1),
taking into account (3.2) and the equation $U_1 = U_e = U_a$, we
obtain:

$$U_0 \approx U_2 \\ \approx U_1 \cdot \frac{R_2}{R_1} = U_a \cdot \frac{R_2}{R_1}$$

Important:

Select resistors so that the voltage drop at R_1 does not exceed
10 V.

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$$I = \frac{U_a}{R}$$

Principle (see fig. 2.3):

By analogy with the ohmic voltage divider (para 3.1.2), the

following applies, taking $C = \frac{1}{\frac{1}{C_1} + \frac{1}{C_2}}$

into account:

$$U_1 = \frac{C_2}{C_1} \cdot U_2$$

and

$$U_2 = \frac{C_1}{C_2} \cdot U_1$$

When $C_1 \gg C_2$ we obtain, for the total voltage

$$U_0 = U_1 + U_2 \quad \text{with } U_1 \ll U_2$$

$$U_0 = U_2$$

$$U_0 = U_1 \cdot \frac{C_1}{C_2}$$

$$= U_a \cdot \frac{C_1}{C_2}$$

Important:

Select capacitors so that the voltage drop at C_1 does not exceed 10 V; for C_2 use capacitors with suitable dielectric strength only. Observe the value printed on the element!

STE capacitors required:

100 pF, 630 V	(578 22)
1 nF, 630 V	(578 25)
10 nF, 100 V	(578 10)
100 nF, 100 V	(578 31)
1 µF, 100 V	(578 15)

3.3 Current measurement

Principle (see fig. 3):

If a current I flows through resistor R a voltage drop U_e occurs at R . According to Ohm's law:

$$U_e = R \cdot I$$

The voltage drop U_e at R is therefore a measure of current intensity I :

$$I = \frac{U_e}{R}$$

Since $U_e = U_a$ this relationship also applies for output voltage U_a :

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3.4 Charge measurement

Principle (see fig. 4):

The following relation exists between the charge Q and the voltage U of a capacitor C :

$$Q = C \cdot U_2$$

Since $U_e = U_a$, output voltage U_a is a measurement for charge Q :

$$Q = U_a \cdot C$$

STE capacitors required:

Cat.-No.	Capacity	Q at $U_A = 1\text{ V}$
578 15	1 μF	10^{-6} As
578 31	100 nF	10^{-7} As
578 10	10 nF	10^{-8} As
578 25	1 nF	10^{-9} As

Table 2

For quantitative investigation of charges in electrostatics, use the Faraday cup (546 12), which is connected by plug (e.g. from 501 561); immerse the specimen - e.g. friction rod - in the field-free area in the cup to record total charge according to fig. 6.

3.5 Offset compensation

Compensate for an offset voltage, which may occur in rare cases with a shorted amplifier input ($U_e = 0$), by zeroing the screwdriver potentiometer ④.


Fig. 5

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Arrangement for demonstrating the ionizing effect of combustion gases or α radiation applied to an open ionization chamber consisting of a zinc plate (546 31) and a grid electrode (546 33).


Fig. 6

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Arrangement for charge measurement in electrostatics.
Important: connect the experimenter to equipment earth (reference potential) via connecting rod (532 16)!