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Using these parts, you can carry out qualitative electrostatic experiments (attraction/repulsion) in conjunction with the Faraday's cup (e.g. 546 12), in which the charge created by rubbing the screen is transferred, and a charge meter (e.g. field electrometer, 540 10, or electrometer amplifier, 532 14, with STE-capacitor, 578 31).

### 1 Description

The glass tubes and the plastics rods are all about 25 cm long. The glass tube (541 02) and plastic rod (541 05) have a small "pan" worked into them in the middle of their length, which is used to fit the tube on the base for the magnetic needle (513 51). The leather (541 21) is soft chamois leather.

### 2 Use

In order to show the repulsion of like and the attraction of unlike charges, it is necessary to provide two of the four tubes and rods with electric charges at a time. For this purpose the tubes and rods must be rubbed. The glass tubes should be rubbed with the leather. Positive charges remain on the glass tube. Instead of the leather dry newsprint can be used. The plastic rods must be rubbed with wool, in order to generate negative charges on the plastic rod.

First, the glass tube with bearing pan is rubbed, and then placed on the needle bearing of the base (513 51). If a rubbed glass tube is now brought into proximity of the first tube, repulsion can be observed; with a rubbed plastic rod, this process demonstrates attraction.

For these experiments the tubes and rods must be clean and dry. If necessary, they should be cleaned with alcohol and a lint-free cloth and dried in hot air (do not clean the plastic rods with fluids containing solvents!).