Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

9/93-Kem/Sf-


Fig. 1

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Neon Franck-Hertz Tube
Holder with Socket and Screen

The neon Franck-Hertz tube, which is operated at room temperature, is used to investigate the discrete energy absorption of neon atoms upon collision with free electrons by observing the electron emission current (compare electron collision experiment by Franck and Hertz using mercury atoms, which can be carried out with the Franck-Hertz tube 555 85). The initial excitation states of Neon are in the range from 16 up to 18eV above the fundamental state. Transitions between these states can be observed as corresponding light bands.

1 Safety instructions


⚠ When using the regulated power supply (52235) alterations may only be made to the circuit with the device switched off, as this device produces contact hazardous voltage. Use safety connecting leads.

⚠ The danger of implosion exists for the tube when it is dropped or banged against other objects, etc.! (internal pressure approx. 10 hPa)

When light phenomena appear outside the region between the grids, reduce the acceleration voltage U_2 and the control voltage U_1 .

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
- ① Cathode heater
- ② Cathode
- ③ Protective cylinder
- ④ Control grid
- ⑤ Acceleration grid
- ⑥ Target electrode

Low-pressure tube with electrode system arranged parallel
 Neon filling: approx. 10 hPa

Heating voltage U_f : 6.3 V (approx. 260 mA)
 Control voltage U_1 between
 cathode and control grid: 5 V
 Acceleration voltage U_2 between
 control grid and acceleration grid: 0...80 V
 Counter voltage U_3 between
 acceleration grid and
 target electrode: 0...10 V

Operation at room temperature
 Excitation level of the Ne-atoms: 16 – 18 eV

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- ① Tube socket
- ② Protective spring
- ③ Connection for the collector

- ④ Principle circuit diagram

Connection for assembly with the Franck-Hertz operating device (555 88):

⑤ DIN socket for the connection cable (555 872)

Connections for assembly with individual devices:

- ⑥ 4-mm sockets for heating voltage U_f
- ⑦ 4-mm sockets for control voltage U_1
- ⑧ 4-mm sockets for accelerating voltage U_2
- ⑨ 4-mm sockets for counter voltage U_3
- ⑩ BNC socket for connection of a current-sensitive measurement configuration for the collector current I

Table-top device, which can also be mounted in an experiment frame.

Dimensions: 20 cm x 29.7 cm x 17 cm

3 Operation

3.1 Operating and measuring equipment

Assembly with individual components:

Power supplies:

- 1 Regulated power supply, 0 up to 300 V .522 35
- 1 Rotary potentiometer, 100 k Ω 537 85 (recommended)
- 1 DC power supply unit, 0 up to ± 15 V .522 30

Measuring amplifier for 10^{-9} A:

- 1 Electrometer amplifier532 14
- 1 Voltage supply 12 V AC, e.g.562 791
- 1 STE resistor 1 G Ω 577 02
- 1 Screened cable BNC/4 mm575 24
- or
- 1 I-measuring amplifier D532 00
- 1 Measurement cable BNC/BNC501 02

Display instruments:

- 2 Voltage measuring instruments (MR 100 V and 10 V) e.g. 531 94 or
- 1 XY-recorder, e.g.575 662
- or (for computer-supported experiment with a DOS-computer)
- 1 CASSYpack E524 007
- 1 Differential amplifier box524 039
- 1 Software "Measuring and Evaluating" 524 112
- 1 Measurement resistor 100 k Ω 536 25
- 1 Measurement resistor 10 k Ω 536 19
- 1 STE capacitor 1 μ F578 35 (recommended)

Connecting leads (follow safety instructions!)

Set-up using the Franck-Hertz operating device:

- 1 Franck-Hertz operating device555 88
- 1 Connection cable555 872

Note: an additional measuring device is not absolutely necessary for the point-by-point recording of the curve described in the practical experiment

Display equipment:

- 2 voltage measuring instruments, MR 100 V and 10 V, e.g.531 94
- or
- 1 XY-recorder, e.g.575 662

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or (for computer-supported experimenting with DOS computers)

- 1 CASSYpack E 524 007
- 1 Software "Measuring and Evaluating" 524 112
- or
- 1 Oscilloscope 575 21
- 2 Screened cables BNC/4 mm 575 24
- Connecting leads

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Fig. 4
Set-up for point-by-point recording of the neon Franck-Hertz curve; for recording the curve connect the XY-recorder or CASSY with the voltage divider circuit instead of the measuring instruments.

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Fig.5
Set-up for investigating the neon Franck-Hertz curve on the oscilloscope using the operating device (555 88).

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
3.3 Experiment procedure

Set-up the arrangement according to Fig. 4 or Fig. 5, but without connecting the recorder yet.

Apply voltages exactly as stated in the described sequence:

- 1.) Set the voltage connected at the rotary potentiometer to 80 V.
- 2.) Switch on the heating current and let the cathode warm up for approx. 1 minute.
- 3.) Set the counter voltage U_3 to approx. 7 V.
- 4.) Set the acceleration voltage U_2 to approx. 70 V.
- 5.) Slowly increase the control voltage U_1 . In the process observe the space between the two grids.
Important: The control voltage which is correct for the tube and the current cathode temperature is achieved, when three (red coloured) light bands can be observed.
- 6.) Set acceleration voltage U_2 to 0 V and then slowly increase it up to the 1st minimum of the collector current (U_2 approx. 20 V).
- 7.) Set the counter voltage U_3 so that collector current almost disappears.

The arrangement is now correctly set for the experiment procedure.

If necessary connect the recorder:
Calibration: X-axis: 10 V/cm var.
Y-axis: 1 V/cm var.

or
Connect CASSY with dU-box via the voltage divider circuit:

Settings: U_2 at input A with dU-box; 0...10 V
 I at input C; calibration 1 V = 1 nA; 0...3 nA
(or 10 nA)

For the experiment procedure slowly increase the acceleration voltage U_2 and observe the dependency of the collector current I .

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