11/95-kem-

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ELM Basic Set ELM Supplementary Set Basic Machine Unit

The basic set (563 480) of the ELM Electrical Machine Teaching Model system permits experiments on electric motors and generators covering a wide variety of configurations with permanent magnet and electromagnet stators and rotors. Depending on the type, the machines are operated either with DC or AC current.

The supplementary set (563 481) additionally makes possible experiments on self-starting motors, with special emphasis on three-phase machines.

All configurations realized with the components of the equipment sets, as well as with (other) individual units, are assembled on the basic machine unit (727 81).

Literature:

Experiments with the electrical machine teaching models (563 492)

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1 Safety notes

- Do not touch rotating parts of the motor or any parts of the experiment setup which the motor drives. Make sure that long hair, scarves, etc. cannot become caught in the rotating parts.
- Apply the supply voltage for DC and AC motors only via the supply sockets ① (Fig. 3) of the basic machine unit.
- Maximum permissible supply voltage: $U_{\max} = 10 \text{ V}$
- Maximum permissible current for carbon brushes:
 $I_{\max} = 1.5 \text{ A}$
- Lubricate the axle frequently - if necessary, clean first with a lint-free cloth - in order to prevent the rotors from seizing up.
- When attaching pole pieces, always use the centering disk (563 17, Fig. 1 ⑧) to prevent subsequent damage from the rotors.
- In assemblies with the two-pole rotor, avoid short-circuits via the brushes (Fig. 7).
- In ELM experiment setups with connected inductances (e.g. DC motor), the limit values for emission of conducted interference in low-voltage networks set forth in European standards and guidelines may easily be exceeded in long-term operation.

The interference caused by the voltage peaks occurring at the inductances (coil rotors with collectors) are initially limited by built-in bipolar suppressor diodes (22 V) connected in parallel. Additional high-frequency energy components are absorbed by a filter at the supply sockets ① provided for motor assemblies (see Fig. 3) on the basic machine unit. In addition, the supply voltage should be limited to 10 V.

These measures are solely for compliance with legal requirements. They do not influence the experiment procedure or results in any way.

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
Fig. 1

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Basic machine unit (727 81)

- ① 2 pairs of 4-mm supply sockets with filter for connecting the voltage for motor experiments;
Filter, consisting of two bar-core chokes $56 \mu\text{H}$ with $10 \text{ k}\Omega$ resistors connected in parallel, for absorbing high-frequency energy components (cf. safety notes) caused by sparking at the collectors.
- ② 3 pairs of 4-mm connection sockets
- ③ Rotor axle: $100 \text{ mm} \times 8 \text{ mm dia.}$
Dimensions of base plate: $200 \text{ mm} \times 297 \text{ mm}$
Weight: 2.5 kg

2.2 Technical data

Coil, 250 turns (563 11)

- Resistance R : approx. 1.5Ω
- Impedance Z : approx. 7Ω
(with coil pole piece on the basic machine unit, $f = 50 \text{ Hz}$)
- Max. current I_{max} : 1.5 A

Two-pole rotor (563 22)

- Number of turns: 2×380
- Resistance R : approx. 1.5Ω
- Impedance $Z (f = 50 \text{ Hz})$: approx. 5.5Ω
- Max. current I_{max} : 1.5 A
- Commutator: two-pole
- Coil termination: 1 suppressor diode 22 V for suppressing voltage peaks (cf. safety notes)

Three-pole rotor (563 23)

- Number of turns: 3×340
- Resistance R : approx. 1.5Ω
- Impedance $Z (f = 50 \text{ Hz})$: approx. 7Ω
- Max. current I_{max} : 1.5 A
- Commutator: three-pole
- Coil termination: 3 suppressor diodes 22 V for suppressing voltage peaks (cf. safety notes)

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Attaching the pole pieces with magnets (identical to procedure for coil pole pieces)


Fig. 5

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Securing the rotor with the brush holder

3 Operation

3.1 Additionally required

Setup aids:

| | |
|--------------------------------|---------|
| Demonstration-experiment frame | 301 300 |
| or | |
| CPS-panel frame C 50 | 666 425 |

Generator drive:

| | |
|--------------------------------|------------------|
| Hand-cranked gear | 563 301 |
| or | |
| Driving unit with driving belt | 727 88 685 96 |

Power supply unit:

| | |
|-----------------------------------|--------|
| Variable low-voltage transformer | 521 39 |
| Three-phase extra-low transformer | 521 29 |

Measuring instrument:

| | |
|---|---------|
| e.g. AV-meter | 531 94 |
| or | |
| Oscilloscope, e.g. | 575 211 |
| or | |
| for computer-assisted recording of measured values (with MS-DOS computer), e.g. | |
| CASSYpack E | 524 007 |
| and  "Measuring and Evaluating" | 524 112 |
| or  "Universal Data Acquisition" | 525 032 |

3.2 Setup

- Place the centering disk on the rotor axle.
- Attach the disk magnets with different poles down on the bare-metal side of a pole piece with magnet and insert the screws in the holes as shown in Fig. 4.
- Screw the pole piece on the centering disk in the desired position on the basic machine unit using the Allan wrench (Fig. 5).
- Proceed in the same way with the coils and the coil pole pieces.
- Remove the centering disk and oil the rotor axle.
- Place the rotor on the rotor axle. Turn the rotor by hand to check whether it runs easily; if necessary, put the centering disk back in the assembly and correct the position of the pole pieces.
- Screw on the brush holder carefully to secure the rotor (also the magnet rotor!) (Fig. 6).

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- When using the commutator rings of the two-pole rotor, align the brush holder parallel to the magnetic field of the stator (Fig. 6).
- Attention: In motor assemblies with two-pole rotor, turn the rotor so that the carbon of the brushes no longer bridges the interrupter gaps of the commutator before connecting the voltage to the brushes (Fig. 7), so that no short circuit can occur.
- To wind a cord on the belt pulley - e.g. for lifting a weight - pull the cord through the loop as shown in Fig. 9 and turn the rotor in the direction of the arrow.


Fig. 7

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Attention: short-circuit when two-pole rotor is at rest. Turn the two-pole rotor before applying the voltage!


Fig. 8

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Brush connection using the slip rings of the coil rotors


Fig. 9

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DC motor

Important: connect the voltage only via sockets ① (see also Fig. 3 and section 1).

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Note to arrangements shown in Fig. 11/12:
Use a CPS-lank panel (666 467) as a spacer.


Fig. 11
Generator mit Handantrieb
Manually driven generator (with hand-cranked gear)


Fig. 12
Drehstromgenerator (Außenpolmaschine)
Three-phase generator (stationary-field machine)