Abb. 1

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
Abb. 3

Fig. 1

The ring-shaped discharge tube is used to demonstrate electrodeless discharge in a He/Ne-filled glass sphere placed in a ring-shaped coil supplied by an alternating current of high frequency. It enables to verify the closed electric flux lines surrounding the magnetic field of a circular coil.

1 Description

The standard equipment includes:

- 1 Discharge tube: glass sphere, approx. 10 cm dia., filled with a He/Ne gas mixture
- 1 Ring-shaped coil: 4 turns, on plastic support into which the glass sphere is fitted

Dimensions:

Base area of the support: approx. 14 cm x 10 cm

Total height: approx. 13 cm

2 Experimental Set-up

2.1 Assembly with high-frequency generator

Apparatus:

Cat. No.

- | | |
|--|-------------------|
| 1 Ring-shaped discharge tube on support | 587 41 |
| 1 High-frequency generator | 587 40*) |
| 1 Power supply unit for tube experiments | 522 35 |
| 1 Plastic rod | 541 04 |
| Connecting leads, 25 cm long | from 501 20 to 24 |

Arrange the assembly according to Fig. 2.

*) No longer available.

Ring-Shaped Discharge Tube on Support


Fig. 2

Connect the high-frequency generator via 6-core cable to the power supply unit and the ring-shaped coil to the high-frequency generator output.

Set selector switch ① on power supply unit (522 35) to lower position and switch on the power supply unit. Set high frequency selector ② on high-frequency generator to a range of approx. 4 to 10 MHz (continuous frequency adjustment via ③).

Ignition of the gas discharge is facilitated by approaching the tube by a rubbed plastic rod.

2.2 Set-up with subdivided spark gap

<i>Apparatus:</i>	<i>Cat. No.</i>
1 Ring-shaped discharge tube on support.....	587 41
1 Leyden jar	544 02
1 U-core with yoke	562 11
1 Clamping device	562 12
1 Mains coil, 500 (or 250*) turns	562 21
	(or 562 22*)
1 Coil, 23 000 turns	562 17
1 Subdivided spark gap	587 46
1 Rheostat, 320 Ω (or 110 Ω*)	537 23
	(or 537 24*)
1 Measuring junction box	502 05
1 Cut-out switch	504 45
Connecting leads	from 501 20ff

*) for 110 V a. c.

3 Explaining the build-up of the discharge

The high-frequency magnetic field inside the tube is linked to a curl field, according to Maxwell's 2 nd equation

$$\oint E \cdot ds = - \frac{d\Phi}{dt}$$

or

$$\text{rot } E = \frac{\delta B}{\delta t}$$

whose field strength E increases with increasing distance from the coil centre. Therefore, the gas-discharge tube brought into the coil field is ignited near the edges.

The ring-shaped gas discharge shows the closed course of the electric flux lines.

Set-up the arrangement according to Fig. 3

Caution! Be extremely careful!

During operation the secondary side carries high voltage!

Keep the experimentation time short, as otherwise the subdivided spark gap will be damaged.

After having carefully checked the set-up, switch on the circuit via the cut-out switch.

Vary the value of dropping resistor until a constant, well visible ring-shaped discharge can be observed.


Fig. 3