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All of the experiments in the area of electrical machine teaching models containing rotating rotors are assembled with the aid of this basic machine unit.

### 1 Safety instructions

- Refrain from reaching into rotating parts of the motor and the parts moved by it in the experiment setup. Make sure that long hair, shawls and similar cannot get caught on rotating parts.
- Apply supply voltage for DC and AC motors exclusively via the supply sockets ① (Fig.1) of the basic machine unit.
- Maximum permissible supply voltage:  
 $U_{\max} = 10 \text{ V}$
- Maximum permissible current for the carbon brushes:  
 $I_{\max} = 1.5 \text{ A}$
- Frequently oil the shafts – required after previous cleaning with a non-fraying rag –, to prevent jamming of the rotors.
- Always use the centering disk (563 17, Fig. 1 ③) when attaching the pole pieces to avoid subsequent damage by the rotor.
- Avoid short-circuits via the carbon brushes in assemblies using the two-pole rotor (Fig. 5).
- When conductors (like e.g. DC motors) are connected in the electrical machine teaching model experiment setups the limiting values regarding the generation of line-bound interference signals in the low-voltage mains as set forth by European standards and stipulations can be slightly exceeded when the system is put into continuous operation.

The interference caused by the voltage peaks brought about by the inductors (coil-type rotors with collectors) is initially limited by the built-in interference suppressors. Other high-frequency energy components are absorbed in a filter in the voltage supply sockets ① (s. Fig. 1) designed for motor assemblies and located on the basic machine unit. In addition to this the supply voltage should be limited to a maximum 10 V.

These measures serve exclusively to fulfill the legal stipulations; they have no effect on the experiment procedure or the experiment results obtained.

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- ① 2 pairs of 4-mm power supply sockets with filter to connect voltage for the motor experiments;

Filter consisting of two rod-type core chokes  $56 \mu\text{H}$  with resistors  $10 \text{ k}\Omega$  connected in parallel for the absorption of high-frequency energy components (see safety instructions), which arise from sparking at the collectors.

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### 3.1 Assembly example

- Place the centering disk on the rotor shaft.
- Attach the disk magnets with uneven pole sides respectively on the blank metal side of a magnetic pole piece and insert the screws into the holes as shown in Fig. 2.
- Tighten the pole piece to the centering disk in the desired position on the basic machine unit using the Allen wrench (Fig. 3).

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Fig. 2: Attachment of the magnetic pole piece (also corresponding to the pole piece for coils)

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Fig. 3: After positioning with the centering disk screw the pole pieces on tightly; oil the shaft.

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- Proceed accordingly with the coils and pole pieces for coils.
- Remove centering disk and oil the rotor shaft.
- Place the rotor on the rotor axis. Manually turn the rotor and check whether it rotates unimpeded; if necessary correct position of the pole piece by reinserting the centering disk.
- Carefully tighten the brush holder to secure the rotor (also the magnetic rotor!) (Fig. 4).

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Fig. 5: Attention: Short-circuit when the two-pole rotor is at a standstill. Turn the two-pole rotor before applying the voltage!

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Fig. 6: Brush connection when using the slip-rings of the coil-type rotors

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- When using the commutator ring of the two-pole rotor, the brush holder is aligned parallel to the magnetic field of the stator (Fig. 4).
- Attention! In the case of motor assemblies with the two pole rotor only connect voltage to the brushes when the rotor has been turned so that the carbon of the brushes no longer bridge the interruption points (Fig. 5) and thus avoiding any short-circuit.