


Acoustics  
Propagation of sound

Propagation of sound in the form of fluctuations in pressure  
Sensor-CASSY, pressure sensor S  $\pm 70$  hPa

**Object of the experiment**

1. Investigate the Propagation of sound in the form of fluctuations in air pressure.

**Setup****Preparation of Sensor-CASSY:**

- Connect the Sensor-CASSY module to a USB port of a computer.
- Run the CASSY-Lab software.
- Configure the following settings in the "Sensor input settings" window.
  - Measuring range:  $-0.7 \dots +0.7$  hPa
  - Measurements: **Average values**
  - Use the  $\rightarrow 0 \leftarrow$  button to determine the zero point
- Configure the following settings in the "Measurement parameters" window.
  - **Add new set of measurements**
  - Interval: **50  $\mu$ s**
  - Time for measurement: **20 ms**
  - Trigger: **pA1, 0.100 hPa, rising edge**


**Apparatus**

1 Resonance tuning fork.....	414 42
1 Tambourine .....	414 32
1 Pressure sensor S, $\pm 70$ hPa .....	524 066
1 Sensor-CASSY 2.....	524 013
1 CASSY-Lab 2 .....	524 220
1 Extension cable, 15-pin .....	501 11
1 Funnel, Boro 3.3, 150 mm diam.....	602 671
1 Rubber tube, 8 mm diam., conforming to DIN 128 65 . 667 183	
2 Saddle bases.....	300 11
1 Stand rod, 223 mm, 10 mm diam., M6 .....	688 808
1 Universal clamp, 0...80 mm .....	666 555
1 Laboratory stand II.....	300 76

Additionally required:  
1 PC with Windows XP or higher

**Procedure**

- Start measuring by pressing the F9 key.
- Hold the tambourine about 1 m from the funnel beat it gently and observe the change in pressure on the screen.
- Remove the trigger via the "Measuring parameters" window.
- Set up the tuning fork about 30 cm in front of the funnel, strike it and start measuring by pressing the F9 key.
- Observe the changes in pressure again.

**Observation****Evaluation**

After the tambourine is beaten, its skin initially moves in one direction. Molecules of air behind it get squeezed together. This means for a short time there is a space where there are particularly many air molecules (compression of air).

If the membrane then moves in the opposite direction, then behind it there will be a space with particularly few molecules (rarefaction of air).

Sound propagates due to a sequence of compressions and rarefactions (pressure fluctuations) in the air.

Such propagating pressure fluctuations in air are called sound waves.

When a tuning fork is struck, similar pressure fluctuations occur due to the vibration of the tines.