

## Thermal behaviour of bodies

## Solid bodies

Linear expansion of tubes when the temperature rises -  
Linear expansion apparatus D

## Object of the experiment

1. Investigating the longitudinal expansion of copper, steel and glass tubes when they are warmed

## Setup



## Safety note:

When exchanging the hot expansion tubes, touch them with a piece of cloth.

## Apparatus

1 Longitudinal expansion apparatus D .....	381 341
1 Erlenmeyer flask, Boro 3.3, 500 ml, narrow neck .....	664 251
1 Petri dish, 120 mm .....	602 740
1 Rubber stopper, one 7-mm hole, 28/34 mm Ø.....	667 265
1 Connector, straight, 8 mm diam., 10 pieces .....	665 222ET2
1 Silicone tubing, 7 mm diam., 1 m .....	604 433
1 Stand base, V-shape, large.....	300 01
1 Stand rod, 75 cm, 12 mm diam. ....	300 43
2 Leybold multiclips .....	301 01
1 Universal clamp .....	666 555
1 Stand ring with stem, 100 mm diam. ....	666 573
1 Wire gauze, 120 x 120 mm .....	608 120
1 Butane gas burner .....	666 711
1 Butane cartridge, 190 g, set of 3 .....	666 712ET3

## Carrying out the experiment

- Fill approx. 300 ml of water into the Erlenmeyer flask and close the flask with the stopper.
- Put the copper tube into the expansion apparatus, and connect it to the Erlenmeyer flask via the tubing.
- Fasten the copper tube in the spring clip so that the groove on the tube rests on the blade of the pointer.
- Shift the expansion tube until the pointer points to the zero of the scale.
- Use the cartridge burner to warm the water in the flask until it boils.
- Observe the pointer deflection, and take the value read from the scale down.
- Repeat the procedure with the steel and the glass tube.

## Measuring example

Material	Value read from the scale	*change in length $\Delta l$ in mm
Copper	3.4	0.68
Steel	2.2	0.44
Glass	1.8	0.36

\*The full-scale value of 5 corresponds to a change in the length of the expansion tube of 1 mm.

## Evaluation

The change in the length of tubes when their temperature rises depends on the material.

In the case of a copper tube the change in length is greater than in the case of steel or glass tubes.

## Remark:

From the length changes obtained in this experiment the linear expansion coefficients  $\alpha$  of the materials copper, steel, and glass can be determined:  $\alpha = \frac{\Delta l}{l \cdot \Delta \vartheta}$

$\Delta l$ : change in the length of the expansion tube

$l$ : initial length of the expansion tube

$\Delta \vartheta$ : difference between the temperature of the warmed tube  $\vartheta_2$  ( $\vartheta_2 = 100^\circ\text{C}$ ) and the room temperature  $\vartheta_1$

Table values for the materials used in the experiment:

$$\text{Copper: } \alpha = 17 \cdot 10^{-6} \frac{1}{\text{K}}$$

$$\text{Steel: } \alpha = 11 \cdot 10^{-6} \frac{1}{\text{K}}$$

$$\text{Glass: } \alpha = 9 \cdot 10^{-6} \frac{1}{\text{K}}$$