


Physik

Chemie · Biologie

Technik

LEYBOLD DIDACTIC GMBH

4/1990 -Sf-

Gebrauchsanweisung
Instruction Sheet
Instrucciones de Servicio

468 74
469 79

Infrarot-Filter Ultraviolett-Filter

Infrared Filter Ultraviolet Filter

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The filters which absorb visible light are used for experiments for which radiation in the infrared or ultraviolet region of the spectrum is required.

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1 Technical Data and Description

Infrared filter (468 74):

- transparent in the infrared region for wavelengths λ between approx. 930 nm and 3200 nm
- non-transparent for visible light

Ultraviolet filter (469 79):

- transparent for wavelengths between approx. 330 nm and 390 nm with a maximum in the region of the intensive mercury line 366 nm,
- non-transparent for visible light, except a small portion in the dark-red region (760 nm).

Dimensions: 50 mm x 50 mm x 1.8 mm

Characteristic features for distinguishing the filters:

- Infrared filter: edges unground; appearing black in counterlight.
Ultraviolet filter: edges ground, appearing dark-violet in counterlight.

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2 Notes on Experiments

It is advisable to clamp the filter into the holder with spring clips (460 22).

Recommended light source for experiments with UV light:

High pressure mercury lamp (451 15) with E 27 holder (451 19), with multiple

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plug for connection to the universal choke (451 30).

Demonstration of UV- and IR radiation:

Excitation of fluorescent or phosphorescent materials to glow or to continue to glow during or after the exciting radiation with UV light; quenching the phosphorescence by IR light:

Ultraviolet phosphorescent screen (469 42), one half covered with fluorescent luminous paint;

- Zinc sulphide screen (468 72) exhibiting phosphorescence;

Card with 10 luminous paints (469 82) exhibiting fluorescence and phosphorescence.

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